

Air

Part No.	Description	Pack	Price in Rs.
Dusters			
1464959	Fisherbrand* Duster, Air Dust'it; Powerful nonconductive plastic bulb duster eliminates dust, lint, and dirt; cleans optical lenses, slides, and other delicate, difficult-to-clean, inaccessible equipment; 2-1/4 dia. x 5 in.H (5.7 x 12.7cm)	Pack of 3	2674
FBA18428	Fisherbrand* MicroDuster III Valve Assembly, With extension hose 3 ft.; Reusable; Surgical grade	Each	1851
23022523	Fisherbrand* Super Friendly Air'lt Duster, Aerosol; 11 oz.; Blows lint, dust and loose dirt away instantly; gas is triple-filtered to 0.3microm; noncorrosive, nonflammable, dry, odorless, and completely volatile; 5 in. polyethylene snorkel	Case of 6	3682

Aluminium Foil

Part No.	Description	Pack	Price in Rs.
Film and Foil Wrapping			
FBA01213	Fisherbrand* Foil, Aluminium; 0.00045 x 3 in. (0.001 x 7.6cm); 1/2 lb. Roll; Width: 7.6cm; Length: 96m;Thickness: 0.011mm	Each	1602
012131	Fisherbrand* Foil, Aluminium; 0.001 x 6 x 63 in. (0.002 x 15.2 x 160cm); 5 lb., 1000 sheets; Width:15cm;Length: 15cm; Thickness: 0.025mm	Pack of 1000	12442
01213100	Fisherbrand* Foil, Aluminium; Standard-gauge roll; 12 in. W x 25 ft. L; 0.018mm Thick; Individual packaging; Dispenser box with cutter bar	Each	108
01213101	Fisherbrand* Foil, Aluminium; Standard-gauge roll; 12 in. W x 200 ft. L; 0.018mm Thick; Individual packaging; Dispenser box with cutter bar	Each	715
01213102	Fisherbrand* Foil, Aluminium; Standard-gauge roll; 12 in. W x 500 ft. L; 0.018mm Thick; Individual packaging; Dispenser box with cutter bar	Each	1757
01213103	Fisherbrand* Foil, Aluminium; Heavy-gauge roll; 18 in. W x 25 ft. L; 0.024mm Thick; Individual packaging; Dispenser box with cutter bar	Each	199
01213104	Fisherbrand* Foil, Aluminium; Heavy-gauge roll; 18 in. W x 50 ft. L; 0.024mm Thick; Individual packaging; Dispenser box with cutter bar	Each	368
01213105	Fisherbrand* Foil, Aluminium; Heavy-gauge roll; 18 in. W x 500 ft. L; 0.024mm Thick; Individual packaging; Dispenser box with cutter bar	Each	3436

Ampoules

Part No.	Description	Pack	Price in Rs.
Ampoule Accessories			
012161	Fisherbrand* Safesnap* Ampoules, accessories; Ampoule opener; Disposable; Polyethylene; For Ampoule size: 1 to 4mL	Pack of 1000	5377
012162	Fisherbrand* Safesnap* Ampoules, accessories; Ampoule opener; Disposable; Polyethylene; For Ampoule size: 5 to 10mL	Pack of 1000	5409

Anemometers

Part No.	Description	Pack	Price in Rs.
Anemometers and Flowmeters			
FBA01239	Fisherbrand* Anemometers Traceable; Digital anemometer with memory; Includes: Two extension wands and two AA alkaline batteries	Each	101235
FBA01241	Fisherbrand* Anemometers Traceable; Anemometer/Thermometer; Digital; Economy; Measures outside wind speed, air velocity in hoods, gas movement in ducts; Includes: 9V alkaline battery, NIST-traceable certificate, instructions	Each	20220
0666428	Fisherbrand* Anemometer/Thermometer Traceable; Anemometer accuracy: +/-3% of full scale; Thermometer accuracy: +/-1 deg. C; Temperature range: 10.0 to 50.0 deg. C (14.0 to 122 deg. F); Temperature resolution: 0.1 deg. (for both C and F); 7 1/4 x 3 x 1 3/4	Each	13964
0666273	Fisherbrand* Anemometer/Thermometer; Hot Wire; Traceable; With data acquisition software for Windows Lite Version, Traceable to NIS; Ideal for cleanrooms and environmental tests; Miniature glass bead thermistor for greater precision	Each	36583
024015	Fisher Scientific* Meter, Multi-purpose Environmental; Enviro-Meter; Traceable; Functions as light meters, anemometer, humidity meter, and thermometer--all in one unit; Bright LCD; MAX/MIN button; Includes Type-K beaded external thermometer probe	Each	12421
Anemometer Accessories			
0666272	Fisherbrand* Anemometers, accessories; Traceable hot wire anemometer/thermometer DAS-4 data logger accessory; Stores up to 8000 readings from all meters; Readings can be taken at intervals from 1 sec. - 99 hrs.; Reading can be d/I to any PC	Each	29887