

Part No.	Description	Pack	Price in Rs.
19068430	Fisherbrand* Fine-Mesh Nylon Hairnets Diameter: 21 in.	Pack of 100	3364
19070625	Fisherbrand* Fine-Mesh Nylon Hairnets Diameter: 18 in.	Pack of 100	2914
19098910	Fisherbrand* Cap, Bouffant; Latex-Free; Spunbonded Polypropylene; Lightweight; non-linting; feature synthetic band stitch with soft floss for secure fit; Flat packed Double polybag; Dia.: 18 in. (46cm); White	Pack of 100	266
19098912	Fisherbrand* Cap, Bouffant; Latex-Free; Spunbonded Polypropylene; Lightweight; non-linting; feature synthetic band stitch with soft floss for secure fit; Flat packed Double polybag; Dia.: 28 in. (71cm); White	Pack of 100	400
19098913	Fisherbrand* Cap, Bouffant; Latex-Free; Spunbonded Polypropylene; Lightweight; non-linting; feature synthetic band stitch with soft floss for secure fit; Flat packed Double polybag; Dia.: 21 in.; White	Pack of 100	280
19098914	Fisherbrand* Cap, Bouffant; Latex-Free; Spunbonded Polypropylene; Lightweight; non-linting; feature synthetic band stitch with soft floss for secure fit; Flat packed Double polybag; Dia.: 24 in.; White	Pack of 100	304
19098915	Fisherbrand* Cap, Bouffant; Latex-Free; Spunbonded Polypropylene; Lightweight; non-linting; feature synthetic band stitch with soft floss for secure fit; Flat packed Double polybag; Dia.: 21 in. (53cm); Blue	Pack of 100	280
19098916	Fisherbrand* Cap, Bouffant; Latex-Free; Spunbonded Polypropylene; Lightweight; non-linting; feature synthetic band stitch with soft floss for secure fit; Flat packed Double polybag; Dia.: 24 in. (61cm); Blue	Pack of 100	304
Facemasks			
0136192	Fisherbrand* Facemask, With Polyfiltration Medium; Three-ply; Bacterial filtration exceeds 92% at 1.0 micron	Case of 400	2720
18960A	Fisherbrand* Facemasks; With Earloops and Clip; Thin middle layer of meltblown polypropylene filter media lies between spunbonded polypropylene layers ultrasonically welded to reduce particle generation; Filtration efficiency is 99%	Case of 600	5481
18960B	Fisherbrand* Facemasks, With Earloops; Thin middle layer of meltblown polypropylene filter media lies between spunbonded polypropylene inner and blue layers ultrasonically welded to reduce particle generation; Filtration efficiency is 99%	Case of 500	3320
18960C	Fisherbrand* Facemasks; With Ties; Thin middle layer of meltblown polypropylene filter media lies between spunbonded polypropylene inner and outer layers ultrasonically welded to reduce particle generation; Filtration efficiency is 99%	Case of 300	2120

Arsenic Testing

Part No.	Description	Pack	Price in Rs.
FBA01405	Fisherbrand* Testing, Chemical, Arsenic; Arsine Generator Complete Unit	Each	9110
0140510	Fisherbrand* Arsine Generator Replacement Parts Absorber Unit	Each	3855
014052	Fisherbrand* Arsine Generator Replacement Parts Flask; Low extractable; Material: Borosilicate Glass; 24/40 Joints; Capacity: 125mL	Each	3204
014055	Fisherbrand* Arsine Generator Replacement Parts Scrubber Unit; For Arsine Generator	Each	4333
S66625	Fisherbrand* Silver Diethyldithiocarbamate, Powder 25 GR, CAS: (1470-61-7)	Each	23680

Bags

Part No.	Description	Pack	Price in Rs.
Autoclaving Bags			
018158	Fisherbrand* Bags, Autoclaving; Elastic Closures High-tensile-strength; Compatible; Remain intact; 200/Pack	Pack of 200	424
0181250	Fisherbrand* Pouches, Sterilization; Instant sealing; 3-1/2 x 5 in.; Peel off protective strip; triple heat-processed seams; 2mil blue transparent polyester/polypropylene laminate film side	Pack of 200	600
0181251	Fisherbrand* Pouches, Sterilization; Instant sealing; 3-1/2 x 9 in.; Peel off protective strip; triple heat-processed seams; 2mil blue transparent polyester/polypropylene laminate film side	Pack of 200	792
0181254	Fisherbrand* Pouches, Sterilization; Instant sealing; 5-1/2 x 10 in.; Peel off protective strip; triple heat-processed seams; 2mil blue transparent polyester/polypropylene laminate film side	Pack of 200	1216
0181255	Fisherbrand* Pouch, Instant sealing sterilization; 7 x 13 in.; Peel off protective strip; Triple heat-processed seams; 2mL-thick blue transparent polyester/polypropylene laminate film side	Pack of 200	2072
0181257	Fisherbrand* Pouches, Sterilization; Instant sealing; 10 x 15 in.; Peel off protective strip; triple heat-processed seams; 2mil blue transparent polyester/polypropylene laminate film side	Pack of 100	1960
0181258	Fisherbrand* Pouches, Sterilization; Instant sealing; 12 x 18 in.; Peel off protective strip; triple heat-processed seams; 2mil blue transparent polyester/polypropylene laminate film side	Pack of 100	2232