

Part No.	Description	Pack	Price in Rs.
19068430	Fisherbrand* Fine-Mesh Nylon Hairnets Diameter: 21 in.	Pack of 100	3364
19070625	Fisherbrand* Fine-Mesh Nylon Hairnets Diameter: 18 in.	Pack of 100	2914
19098910	Fisherbrand* Cap, Bouffant; Latex-Free; Spunbonded Polypropylene; Lightweight; non-linting; feature synthetic band stitch with soft floss for secure fit; Flat packed Double polybag; Dia.: 18 in. (46cm); White	Pack of 100	266
19098912	Fisherbrand* Cap, Bouffant; Latex-Free; Spunbonded Polypropylene; Lightweight; non-linting; feature synthetic band stitch with soft floss for secure fit; Flat packed Double polybag; Dia.: 28 in. (71cm); White	Pack of 100	400
19098913	Fisherbrand* Cap, Bouffant; Latex-Free; Spunbonded Polypropylene; Lightweight; non-linting; feature synthetic band stitch with soft floss for secure fit; Flat packed Double polybag; Dia.: 21 in.; White	Pack of 100	280
19098914	Fisherbrand* Cap, Bouffant; Latex-Free; Spunbonded Polypropylene; Lightweight; non-linting; feature synthetic band stitch with soft floss for secure fit; Flat packed Double polybag; Dia.: 24 in.; White	Pack of 100	304
19098915	Fisherbrand* Cap, Bouffant; Latex-Free; Spunbonded Polypropylene; Lightweight; non-linting; feature synthetic band stitch with soft floss for secure fit; Flat packed Double polybag; Dia.: 21 in. (53cm); Blue	Pack of 100	280
19098916	Fisherbrand* Cap, Bouffant; Latex-Free; Spunbonded Polypropylene; Lightweight; non-linting; feature synthetic band stitch with soft floss for secure fit; Flat packed Double polybag; Dia.: 24 in. (61cm); Blue	Pack of 100	304
Facemasks			
0136192	Fisherbrand* Facemask, With Polyfiltration Medium; Three-ply; Bacterial filtration exceeds 92% at 1.0 micron	Case of 400	2720
18960A	Fisherbrand* Facemasks; With Earloops and Clip; Thin middle layer of meltblown polypropylene filter media lies between spunbonded polypropylene layers ultrasonically welded to reduce particle generation; Filtration efficiency is 99%	Case of 600	5481
18960B	Fisherbrand* Facemasks, With Earloops; Thin middle layer of meltblown polypropylene filter media lies between spunbonded polypropylene inner and blue layers ultrasonically welded to reduce particle generation; Filtration efficiency is 99%	Case of 500	3320
18960C	Fisherbrand* Facemasks; With Ties; Thin middle layer of meltblown polypropylene filter media lies between spunbonded polypropylene inner and outer layers ultrasonically welded to reduce particle generation; Filtration efficiency is 99%	Case of 300	2120

Arsenic Testing

Part No.	Description	Pack	Price in Rs.
FBA01405	Fisherbrand* Testing, Chemical, Arsenic; Arsine Generator Complete Unit	Each	9110
0140510	Fisherbrand* Arsine Generator Replacement Parts Absorber Unit	Each	3855
014052	Fisherbrand* Arsine Generator Replacement Parts Flask; Low extractable; Material: Borosilicate Glass; 24/40 Joints; Capacity: 125mL	Each	3204
014055	Fisherbrand* Arsine Generator Replacement Parts Scrubber Unit; For Arsine Generator	Each	4333
S66625	Fisherbrand* Silver Diethyldithiocarbamate, Powder 25 GR, CAS: (1470-61-7)	Each	23680

Bags

Part No.	Description	Pack	Price in Rs.
Autoclaving Bags			
018158	Fisherbrand* Bags, Autoclaving; Elastic Closures High-tensile-strength; Compatible; Remain intact; 200/Pack	Pack of 200	424
0181250	Fisherbrand* Pouches, Sterilization; Instant sealing; 3-1/2 x 5 in.; Peel off protective strip; triple heat-processed seams; 2mil blue transparent polyester/polypropylene laminate film side	Pack of 200	600
0181251	Fisherbrand* Pouches, Sterilization; Instant sealing; 3-1/2 x 9 in.; Peel off protective strip; triple heat-processed seams; 2mil blue transparent polyester/polypropylene laminate film side	Pack of 200	792
0181254	Fisherbrand* Pouches, Sterilization; Instant sealing; 5-1/2 x 10 in.; Peel off protective strip; triple heat-processed seams; 2mil blue transparent polyester/polypropylene laminate film side	Pack of 200	1216
0181255	Fisherbrand* Pouch, Instant sealing sterilization; 7 x 13 in.; Peel off protective strip; Triple heat-processed seams; 2mL-thick blue transparent polyester/polypropylene laminate film side	Pack of 200	2072
0181257	Fisherbrand* Pouches, Sterilization; Instant sealing; 10 x 15 in.; Peel off protective strip; triple heat-processed seams; 2mil blue transparent polyester/polypropylene laminate film side	Pack of 100	1960
0181258	Fisherbrand* Pouches, Sterilization; Instant sealing; 12 x 18 in.; Peel off protective strip; triple heat-processed seams; 2mil blue transparent polyester/polypropylene laminate film side	Pack of 100	2232

Part No.	Description	Pack	Price in Rs.
01814A	Fisherbrand* Bag, Autoclave; HDPE; Orange, With sterilization indicator; Printed in English and Spanish w/OSHA-compliant graphics; Max. temp: 125 deg. C; Supplied with wire closure; Size: 14 W x 19 in. L	Pack of 200	2540
01814B	Fisherbrand* Bag, Autoclave; HDPE; Orange, With sterilization indicator; Printed in English and Spanish w/OSHA-compliant graphics; Max. temp: 125 deg. C; Supplied with wire closure; Size: 19 W x 23 in. L	Pack of 200	4241
01814C	Fisherbrand* Bag, Autoclave; HDPE; Orange, With sterilization indicator; Printed in English and Spanish w/OSHA-compliant graphics; Max. temp: 125 deg. C; Supplied with wire closure; Size: 25 W x 35 in. L	Pack of 200	7939
01814D	Fisherbrand* Bags, Autoclaving; HDPE; Orange; With sterilization indicator or closure; Printed in English; Size: 38 x 47 in.	Pack of 100	7854
018151	Fisherbrand* Bag, Autoclave; Printed Orange, High Density, Polyethylene; use at 121 deg. C (250 deg. F), Size: 8-1/2 W x 11 in. L; not supplied with wire tie closures	Pack of 100	414
01815A	Fisherbrand* Bag, Autoclave; Printed Orange, High Density, Polyethylene; use at 121 deg. C (250 deg. F), Size: 12 W x 24 in. L, not supplied with wire tie closures	Case of 100	1048
01815B	Fisherbrand* Bags, Autoclaving; HDPE; Printed in English; Orange; Without sterilization indicator or closure; Size: 24 x 30 in.	Case of 100	2396
01815C	Fisherbrand* Bag, Autoclave; Printed Orange, High Density, Polyethylene; use at 121 deg. C (250 deg. F), Size: 24 W x 36 in. L, not supplied with wire tie closures	Case of 100	2852
018141	Fisherbrand* Bag, Biohazard Autoclave; Polyethylene; Plain; Clear; Max. temp: 125 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); Thickness: 1.5mil; Size: 12 W x 24 in. L	Pack of 200	1883
018142	Fisherbrand* Bag, Biohazard Autoclave; Polyethylene; Plain; Clear; Max. temp: 125 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); Thickness: 1.5mil; Size: 24 W x 30 in. L	Pack of 200	4388
018143	Fisherbrand* Bag, Biohazard Autoclave; Polyethylene; Plain; Clear; Max. temp: 125 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); Thickness: 1.5mil; Size: 24 W x 36 in. L	Pack of 200	5273
01827A	Fisherbrand* Bag, Biohazard Autoclave; Polyethylene; Printed in English and Spanish to conform w/OSHA standard 29CFR1910.1030; Clear; Max. temp: 127 deg. C; Pass Dart Impact Resistance Test; Thickness: 1.5mil; Size: 8 x 12 in.	Pack of 200	887
01827B	Fisherbrand* Bag, Biohazard Autoclave; Polyethylene; Printed in English and Spanish to conform w/OSHA standard 29CFR1910.1030; Clear; Max. temp: 127 deg. C; Pass Dart Impact Resistance Test; Thickness: 1.5mil; Size: 12 x 24 in.	Pack of 200	1856
01827C	Fisherbrand* Bag, Biohazard Autoclave; Polyethylene; Printed in English and Spanish to conform w/OSHA standard 29CFR1910.1030; Clear; Max. temp: 127 deg. C; Pass Dart Impact Resistance Test; Thickness: 1.5mil; Size: 24 x 30 in.	Pack of 200	4399
01827D	Fisherbrand* Bag, Biohazard Autoclave; Polyethylene; Printed in English and Spanish to conform w/OSHA standard 29CFR1910.1030; Clear; Max. temp: 127 deg. C; Pass Dart Impact Resistance Test; Thickness: 1.5mil; Size: 24 x 36 in.	Pack of 200	5230
01829A	Fisherbrand* Bag, Biohazard Autoclave; Polyethylene; Printed in English and Spanish to conform w/OSHA standard 29CFR1910.1030; Red; Max. temp: 127 deg. C; Pass Dart Impact Resistance Test; Thickness: 1.5mil; Size: 8 x 12 in.	Pack of 200	899
01829B	Fisherbrand* Bag, Biohazard Autoclave; Polyethylene; Printed in English and Spanish to conform w/OSHA standard 29CFR1910.1030; Red; Max. temp: 127 deg. C; Pass Dart Impact Resistance Test; Thickness: 1.5mil; Size: 14 x 19 in.	Pack of 200	1781
01829C	Fisherbrand* Bag, Biohazard Autoclave; Polyethylene; Printed in English and Spanish to conform w/OSHA standard 29CFR1910.1030; Red; Max. temp: 127 deg. C; Pass Dart Impact Resistance Test; Thickness: 1.5mil; Size: 19 x 23 in.	Pack of 200	3061
01829D	Fisherbrand* Bag, Biohazard Autoclave; Polyethylene; Printed in English and Spanish to conform w/OSHA standard 29CFR1910.1030; Red; Max. temp: 127 deg. C; Pass Dart Impact Resistance Test; Thickness: 1.5mil; Size: 25 x 35 in.	Pack of 200	5428
01829E	Fisherbrand* Bag, Biohazard Autoclave; Polyethylene; Printed in English and Spanish to conform w/OSHA standard 29CFR1910.1030; Red; Max. temp: 127 deg. C; Pass Dart Impact Resistance Test; Thickness: 1.5mil; Size: 38 x 48 in.	Pack of 100	5278
01829F	Fisherbrand* Bag, Biohazard Autoclave; Polyethylene; Printed in English and Spanish to conform w/OSHA standard 29CFR1910.1030; Red; Max. temp: 127 deg. C; Pass Dart Impact Resistance Test; Thickness: 1.5mil; Size: 33 x 30 in.	Pack of 200	5854
01831A	Fisherbrand* Bag, Biohazard Autoclave; Polyethylene; Printed in English and Spanish to conform w/OSHA standard 29CFR1910.1030; Red; Max. temp: 127 deg. C; Pass Dart Impact Resistance Test; Thickness: 4.0mil; Size: 24 x 36 in.	Case of 100	7451
01831B	Fisherbrand* Bag, Biohazard Autoclave; Polyethylene; Printed in English and Spanish to conform w/OSHA standard 29CFR1910.1030; Red; Max. temp: 127 deg. C; Pass Dart Impact Resistance Test; Thickness: 4.0mil; Size: 38 x 48 in.	Case of 50	7294
01832A	Fisherbrand* Bag, Biohazard Autoclave; Polyethylene; Plain; Clear; Max. temp: 127 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); Thickness: 4.0mil; Size: 24 x 36 in.	Case of 100	7257
22044558	Fisherbrand* Bags, Biohazard; Autoclave; Polyethylene; 1.5mil thick; 12 L x 24 in. W; clear; without indicators	Case of 100	1016
22044559	Fisherbrand* Bags, Biohazard; Autoclave; Polyethylene; Without indicators; Clear; 1.5mil thick; Size: 24 x 30 in.	Case of 100	2269
22044560	Fisherbrand* Bags, Biohazard; Autoclave; Polyethylene; Without indicators; Clear; 1.5mil thick; Size: 24 x 36 in.	Case of 100	2700
22044561	Fisherbrand* Bags, Biohazard; Autoclave; Polyethylene; Without indicators; Clear; 1.25mil thick; Size: 8 x 12 in	Case of 1000	3882

Part No.	Description	Pack	Price in Rs.
22044562	Fisherbrand* Bags, Biohazard; Autoclave; Polyethylene; 1.25mil thick; 14 L x 19 in. W; clear; without indicators	Case of 200	1551
22044563	Fisherbrand* Bags, Biohazard; Autoclave; Polyethylene; Without indicators; Red; 1.25mil thick; Size:19 x 23 in	Case of 200	2813
22044564	Fisherbrand* Bags, Biohazard; Autoclave; Polyethylene; Without indicators; Red; 1.5mil thick; Size: 25 x 35 in	Case of 200	5602
22044565	Fisherbrand* Bags, Biohazard; Autoclave; Polyethylene; Without indicators; Clear; 1.5mil thick; Size: 36 x 45 in.	Case of 100	4921
14375175	Fisherbrand* Holders for Autoclave Bags; High-density polyethylene with lid; Capacity: 8 qt.; Accepts bags with minimum width of 12 in.	Each	1196
14375176	Fisherbrand* Holders for Autoclave Bags; High-density polyethylene with lid; Capacity: 24 qt.; Accepts bags with minimum width of 16 in.	Each	2622
14375177	Fisherbrand* Holders for Autoclave Bags; High-density polyethylene with Lid; Capacity: 48 qt.; Accepts bags with minimum width of 20 in.	Each	3917
018261	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Plain; Clear; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); With elastic closure; Thickness: 2.0mil; Size: 8 x 12 in.	Pack of 200	1149
018262	Fisherbrand* Bag, Biohazard, Autoclave; Polypropylene; Clear; With autoclavable elastic steam-venting closures; Pass Dart Impact Resistance Test of 165g; Size: 10 x 10 in. Thickness: 2mil	Pack of 200	1332
018263	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Plain; Clear; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); With elastic closure; Thickness: 2.0mil; Size: 10 x 12 in.	Pack of 200	1363
018264	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Plain; Clear; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); With elastic closure; Thickness: 2.0mil; Size: 10 x 15 in.	Pack of 200	1610
018265	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Plain; Clear; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); With elastic closure; Thickness: 2.0mil; Size: 12 x 24 in.	Pack of 200	2796
018266	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Plain; Clear; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); With elastic closure; Thickness: 2.0mil; Size: 24 x 30 in.	Pack of 200	6364
018267	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Plain; Clear; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); With elastic closure; Thickness: 2.0mil; Size: 24 x 36 in.	Pack of 200	7619
01826A	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Printed in English/Spanish to conform w/OSHA standard 29CFR1910.1030; Clear; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test; With elastic closure; 2.0mil; Size: 12 x 24 in.	Pack of 200	2618
01826B	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Printed in English/Spanish to conform w/OSHA standard 29CFR1910.1030; Clear; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test; With elastic closure; 2.0mil; Size: 24 x 30 in.	Pack of 200	6306
01826C	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Printed in English/Spanish to conform w/OSHA standard 29CFR1910.1030; Clear; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test; With elastic closure; 2.0mil; Size: 24 x 36 in.	Pack of 200	7550
01826D	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Printed in English/Spanish to conform w/OSHA standard 29CFR1910.1030; Clear; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test; With elastic closure; 2.0mil; Size: 31 x 38 in.	Pack of 200	9905
01826E	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Printed in English/Spanish to conform w/OSHA standard 29CFR1910.1030; Clear; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test; With elastic closure; 2.0mil; Size: 37 x 48 in.	Pack of 100	7446
01828A	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Printed in English/Spanish to conform w/OSHA standard 29CFR1910.1030; Red; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test; With elastic closure; 2.0mil; Size: 8 x 12 in.	Pack of 200	1143
01828B	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Printed in English and Spanish to conform w/OSHA standard 29CFR1910.1030; Red; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test; With elastic closure; 2.0mil; Size: 14 x 19 in.	Pack of 200	2408
01828C	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Printed in English and Spanish to conform w/OSHA standard 29CFR1910.1030; Red; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test; With elastic closure; 2.0mil; Size: 19 x 23 in.	Pack of 200	4047
01828D	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Printed in English and Spanish to conform w/OSHA standard 29CFR1910.1030; Red; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test; With elastic closure; 2.0mil; Size: 25 x 35 in.	Pack of 200	7598
01828E	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Printed in English and Spanish to conform w/OSHA standard 29CFR1910.1030; Red; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test; With elastic closure; 2.0mil; Size: 37 x 48 in.	Pack of 100	7467
01830A	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Plain; Red; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); With elastic closure; Thickness: 2.0mil; Size: 8 x 12 in.	Pack of 200	1143
01830B	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Plain; Red; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); With elastic closure; Thickness: 2.0mil; Size: 14 x 19 in.	Pack of 200	2528
01830C	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Plain; Red; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); With elastic closure; Thickness: 2.0mil; Size: 19 x 23 in.	Pack of 200	4116
01830D	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Plain; Red; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); With elastic closure; Thickness: 2.0mil; Size: 25 x 35 in.	Pack of 200	7545
01830E	Fisherbrand* Bag, Biohazard Autoclave; Polypropylene; Plain; Red; Max. temp: 140 deg. C; Pass Dart Impact Resistance Test of 165g (using ASTM D1709-85); With elastic closure; Thickness: 2.0mil; Size: 37 x 48 in.	Pack of 100	7393

Part No.	Description	Pack	Price in Rs.
Biopsy Bags			
15182500	Fisherbrand* Bags, Biopsy; Nylon; Small; Nylon cloth with precise mesh diameter; 1.0 in.L x 1.6 in.W (25mm x 40mm)	Pack of 1000	21603
15182501	Fisherbrand* Bags, Biopsy; Nylon; Medium; Nylon cloth with precise mesh diameter; 1.8 in.L x 3.0 in.W (45mm x 75mm)	Pack of 200	7761
15182502	Fisherbrand* Bags, Biopsy; Nylon; Large; Nylon cloth with precise mesh diameter; 3.0 in.L x 3.7 in.W (75mm x 95mm)	Pack of 500	16322
15182116	Fisherbrand* Biopsy Bag, Nylon, Small; Allow unrestricted fluid movement around tissue; Autoclaveable; Peel-open bag exposes specimen for rapid removal; Reduced reagent carry-over and pressure artifact; No biohazard label; 30 x 50mm	Pack of 500	26138
15182117	Fisherbrand* Biopsy Bag, Nylon, Small; Allow unrestricted fluid movement around tissue; Autoclaveable; Peel-open bag exposes specimen for rapid removal; Reduced reagent carry-over and pressure artifact; No biohazard label; 45 x 75mm	Pack of 1000	34405
Controlled Environment Bags			
19033616	Fisherbrand* Precision Clean II* Bag; Cleanroom; Clear; 1.5 mil; W x L: 36 x 36 in.	Case of 200	18990
19033619	Fisherbrand* Precision Clean II* Bag; Cleanroom; Clear; 2 mil; W x L: 5 x 8 in.	Case of 2000	13023
19033622	Fisherbrand* Precision Clean II* Bag; Cleanroom; Clear; 2 mil; W x L: 8 x 12 in.	Case of 1500	19992
19033625	Fisherbrand* Precision Clean II* Bag; Cleanroom; Clear; 2 mil; W x L: 12 x 18 in.	Case of 1000	27594
19033626	Fisherbrand* Precision Clean II* Bag; Cleanroom; Clear; 2mil; W x L: 18 x 24 in.	Case of 500	29004
19033628	Fisherbrand* Precision Clean II* Bag; Cleanroom; Clear; 4mil; W x L: 24 x 30 in.	Case of 150	25864
19033631	Fisherbrand* Precision Clean II* Bag; Cleanroom; Clear; 6 mil; W x L: 8 x 10 in.	Case of 500	27056
19033632	Fisherbrand* Precision Clean II* Bag; Cleanroom; Clear; 6 mil; W x L: 10 x 12 in.	Case of 500	23222
19033633	Fisherbrand* Precision Clean II* Bag; Cleanroom; Clear; 6mil; W x L: 12 x 15 in.	Case of 500	33934
19033634	Fisherbrand* Precision Clean II* Bags; Clear; Class ISO 5 (Formerly Class 100); Silicon-free; 2mil; 8 x 13 in.	Case of 1500	13904
19033635	Fisherbrand* Precision Clean II* Bag; Cleanroom; Clear; 6mil; W x L: 18 x 24 in.	Case of 200	31327
Sample Bags			
19075388D	Fisherbrand* Bags; Biohazard Specimen; Polyethylene; Printed with biohazard graphics; Zipper closure; FDA-approved PE, Pocket for paperwork; 2mil-thick pouch; 1.5mil-thick pocket; Size: 6 x 10 in.	Case of 1000	7863
19075388E	Fisherbrand* Bags; Biohazard Specimen; Polyethylene; Printed with biohazard graphics; Zipper closure; FDA-approved PE, Pocket for paperwork; 2mil-thick pouch; 1.5mil-thick pocket; Size: 5 x 9 in.	Case of 1000	5919
19075388F	Fisherbrand* Bags; Biohazard Specimen; Polyethylene; Printed with biohazard graphics; Zipper closure; FDA-approved PE, Pocket for paperwork; 2mil-thick pouch; 1.5mil-thick pocket; Size: 10 x 10 in.	Case of 1000	10489
0180000	Fisherbrand* Bags, Specimen; Polyethylene; Three-wall; With document compartment and printed biohazard label; Comply with OSHA and NCCLS; Size: 4 x 6 in.	Case of 1000	1454
0180001	Fisherbrand* Bags, Specimen; Polyethylene; Three-wall; With document compartment and printed biohazard label; Reclosable zipper seal; Comply with OSHA and NCCLS; Size: 6 x 6 in.	Case of 1000	3456
0180002	Fisherbrand* Bags, Specimen; Polyethylene; Three-wall; With document compartment; Comply with OSHA and NCCLS; Size: 6 x 9 in.	Case of 1000	2677
0180003	Fisherbrand* Bags, Specimen; Zippered; Polyethylene; Three-wall; With attached document compartment; Comply with OSHA and NCCLS guidelines; Printed; L x W: 6 x 9 in. (15.2 x 23cm)	Case of 1000	2703
0180004	Fisherbrand* Bags, Specimen; Zippered; Polyethylene; Three-wall; With attached document compartment; Comply with OSHA and NCCLS guidelines; Printed; L x W: 6 x 9 in. (15.2 x 23cm); With requisition pouch and red STAT labeling	Case of 1000	2899
0180005	Fisherbrand* Bags, Specimen; Zippered; Polyethylene; Three-wall; With attached document compartment; Comply with OSHA and NCCLS guidelines; Printed; L x W: 6 x 9 in. (15.2 x 23cm); With 3 x 5 in. pad for absorbing water, urine, or blood	Case of 500	3807
0180006	Fisherbrand* Bags, Specimen; Polyethylene; Three-wall; With document compartment; Reclosable zipper seal; Comply with OSHA and NCCLS; With Scoreline; Size: 6 x 9 in.	Case of 1000	2963
0180007	Fisherbrand* Bags, Specimen; Polyethylene; Three-wall; With document compartment; Comply with OSHA and NCCLS; Size: 8 x 8 in.	Case of 1000	3229
0180008	Fisherbrand* Bags, Specimen; Zippered; Polyethylene; Three-wall; With attached document compartment; Comply with OSHA and NCCLS guidelines; Printed; L x W: 8 x 10 in. (20 x 25cm)	Case of 1000	3732
0180009	Fisherbrand* Bags, Specimen; Zippered; Polyethylene; Three-wall; With attached document compartment; Comply with OSHA and NCCLS guidelines; Printed; L x W: 12 x 15 in. (30 x 38cm)	Case of 1000	7805

Part No.	Description	Pack	Price in Rs.
0180010	Fisherbrand* Bags, Transport; Color-coded zippers; 6 x 9 in. (15 x 23cm); With document pouch; Printed with OSHA-required biohazard labeling Color: Blue	Case of 1000	2899
0180011	Fisherbrand* Bags, Transport; Color-coded zippers; 6 x 9 in. (15 x 23cm); With document pouch; Printed with OSHA-required biohazard labeling Color: Green	Case of 1000	2868
0180012	Fisherbrand* Bags, Transport; Color-coded zippers; Document pouch keeps paperwork separate from specimens; Printed with OSHA-required biohazard labeling; Color: Orange; 6 x 9 in. (15 x 23cm)	Case of 1000	2899
0180013	Fisherbrand* Bags, Transport; Color-coded zippers; 6 x 9 in. (15 x 23cm); With document pouch; Printed with OSHA-required biohazard labeling Color: Red	Case of 1000	2868
0180014	Fisherbrand* Bags, Transport; Color-coded zippers; 6 x 9 in. (15 x 23cm); With document pouch; Printed with OSHA-required biohazard labeling Color: Yellow	Case of 1000	2868
0180015	Fisherbrand* Bags, Transport; Color-coded zippers; 6 x 9 in. (15 x 23cm); With document pouch; Printed with OSHA-required biohazard labeling Color: Assorted	Case of 1000	4673
18028321	Fisherbrand* Bags, Specimen Transport; Infecon; For collecting, storing, transporting tissue samples; Document pouch keeps paperwork with specimen; Printed with Biohazard Symbol; Double-lock closure makes it liquid tight; 8W x 10 in. L	Pack of 250	14132
19028320	Fisherbrand* Bags, Specimen Transport; Infecon; For collecting, storing, transporting tissue samples; Document pouch keeps paperwork with specimen; Printed with Biohazard Symbol; Double-lock closure makes it liquid tight; 12W x 12 in. L	Pack of 250	15627
19075210	Fisherbrand* Bags, Specimen Transport; Infecon; For collecting, storing, transporting tissue samples; Document pouch keeps paperwork with specimen; Printed with Biohazard Symbol; Double-lock closure makes it liquid tight; 16W x 16 in. L	Case of 100	8755
19287215	Fisherbrand* Bags, Specimen Transport; Infecon; For collecting, storing, transporting tissue samples; Document pouch keeps paperwork with specimen; Printed with Biohazard Symbol; Double-lock closure makes it liquid tight; 6W x 9 in. L	Case of 1000	39425
22287216	Fisherbrand* Bags, Specimen Transport; Infecon; For collecting, storing, transporting tissue samples; Document pouch keeps paperwork with specimen; Printed with Biohazard Symbol; Double-lock closure makes it liquid tight; 9W x 12 in. L	Case of 1000	60706
19075206	Fisherbrand* Bags, Specimen Storage Bitran; PE; For short-term storage; For collecting, transporting, storing tissue specimens; For use as mini-anaerobic environments; Resealable double-lock zipper closure; Size: 16 x 16 in.	Case of 100	7785
19075207	Fisherbrand* Bags, Specimen Storage Bitran; PE; For short-term storage; For collecting, transporting, storing tissue specimens; For use as mini-anaerobic environments; Resealable double-lock zipper closure; Size: 24 x 24 in.	Case of 100	16369
19075208	Fisherbrand* Bags, Specimen Storage Bitran; Saranex odor ctrl film; For long-term storage; For collecting, transporting, storing tissue specimens; Use as mini-anaerobic environments; Resealable double-lock zipper closure; Size: 16 x 16 in.	Case of 100	12011
19075209	Fisherbrand* Bags, Specimen Bitran; For collecting, transporting and storing tissue specimens; Resealable double-lock zipper closure; Clear; Saranex; Long-term storage, with odor control; 24 x 24 in.	Case of 100	26829
19240093	Fisherbrand* Bags, Specimen Storage Bitran; Saranex odor ctrl film; For long-term storage; For collecting, transporting, storing tissue specimens; Use as mini-anaerobic environments; Resealable double-lock zipper closure; Size: 2 x 4 in.	Pack of 250	5354
19240101	Fisherbrand* Bags, Specimen Storage Bitran; Saranex odor ctrl film; For long-term storage; For collecting, transporting, storing tissue specimens; Use as mini-anaerobic environments; Resealable double-lock zipper closure; Size: 3 x 6 in.	Pack of 250	7424
19240119	Fisherbrand* Bags, Specimen Bitran; For collecting, transporting and storing tissue specimens; Resealable double-lock zipper closure; Clear; Saranex; Long-term storage, with odor control; 6 x 6 in.	Pack of 250	9272
19240127	Fisherbrand* Bags, Specimen Storage Bitran; Saranex odor ctrl film; For long-term storage; For collecting, transporting, storing tissue specimens; Use as mini-anaerobic environments; Resealable double-lock zipper closure; Size: 7 x 8 in.	Pack of 250	11235
19240143	Fisherbrand* Bags, Specimen Bitran; For collecting, transporting and storing tissue specimens; Resealable double-lock zipper closure; Clear; Saranex; Long-term storage, with odor control; 9 x 12 in.	Pack of 250	16666
19240150	Fisherbrand* Bags, Specimen Storage Bitran; Saranex odor ctrl film; For long-term storage; For collecting, transporting, storing tissue specimens; Use as mini-anaerobic environments; Resealable double-lock zipper closure; Size: 12 x 12 in.	Pack of 250	18629
19240168	Fisherbrand* Bags, Specimen Storage Bitran; PE; For short-term storage; For collecting, transporting, storing tissue specimens; For use as mini-anaerobic environments; Resealable double-lock zipper closure; Size: 3 x 6 in.	Pack of 250	5776
19240192	Fisherbrand* Bags, Specimen Storage Bitran; PE; For short-term storage; For collecting, transporting, storing tissue specimens; For use as mini-anaerobic environments; Resealable double-lock zipper closure; Size: 6 x 6 in.	Pack of 250	7531
19240200	Fisherbrand* Bags, Specimen Storage Bitran; PE; For short-term storage; For collecting, transporting, storing tissue specimens; For use as mini-anaerobic environments; Resealable double-lock zipper closure; Size: 7 x 8 in.	Pack of 250	8427
19240218	Fisherbrand* Bags, Specimen Storage Bitran; PE; For short-term storage; For collecting, transporting, storing tissue specimens; For use as mini-anaerobic environments; Resealable double-lock zipper closure; Size: 9 x 12 in.	Pack of 250	11906
19240226	Fisherbrand* Bags, Specimen Bitran; For collecting, transporting and storing tissue specimens; Resealable double-lock zipper closure; Clear; Polyethylene; Short-term storage; 12 x 12 in.	Pack of 250	14382

Part No.	Description	Pack	Price in Rs.
FBA01823	Fisherbrand* Specimen Bag, Lab-Guard, with Biohazard Label, 6 x 9 in. (15 x 23 cm),	Case of 1000	5169
22306297	Fisherbrand* Lab-Guard Specimen Bags 6x9 in. Plain M/Cs.	Case of 1000	4717
22310302	Fisherbrand* Lab-Guard Specimen Bags 12x15in Biohaz M/Cs.	Case of 1000	16289
1200915A	Fisherbrand* Bags, controlled environment; Bags, sample; Low-density polyethylene; Transparent; Designed to fit Fisherbrand dry waste containers; Capacity: 6.25 gal.; Size: 28L x 22W in. (71 x 56cm)	Case of 100	2232
19075388A	Fisherbrand* Bags, controlled environment; Bags, sample; Flat; Polyethylene; For short-term storage; Thickness: 4mil; Size: 12L x 8 in. W (30.5 x 20.3cm)	Pack of 100	1047
19075388B	Fisherbrand* Bags, controlled environment; Bags, sample; Flat; Polyethylene; For short-term storage; Thickness: 4mil; Size: 10L x 6 in. W (25.4 x 15.2cm)	Pack of 100	811
19075388C	Fisherbrand* Bags, controlled environment; Bags, sample; Flat; Polyethylene; For short-term storage; Thickness: 4mil; Size: 16L x 10 in. W (40.6 x 25.4cm)	Pack of 100	1447
19075388G	Fisherbrand* Bags, controlled environment; Bags, sample; Flat; Polyethylene; For short-term storage; Thickness: 4mil; Size: 6L x 4 in. W (15.3 x 10.2cm)	Pack of 100	501
19075388H	Fisherbrand* Bags, controlled environment; Bags, sample; Flat; Polyethylene; For short-term storage; Thickness: 4mil; Size: 18L x 12 in. W (45.7 x 30.5cm)	Pack of 100	1909
19075388J	Fisherbrand* Bags, controlled environment; Bags, sample; Flat; Polyethylene; For short-term storage; Thickness: 4mil; Size: 12L x 10 in. W (30.5 x 25.4cm)	Pack of 100	1191
19075388K	Fisherbrand* Bags, controlled environment; Bags, sample; Flat; Polyethylene; For short-term storage; Thickness: 4mil; Size: 20L x 12 in. W (50.8 x 30.5cm)	Pack of 100	2059
0181210C	Fisherbrand* Bags, controlled environment; Bags, polyethylene; 4mil thick; 6L x 12 in. W (15.2 x 30.5cm)	Case of 300	1786
0181210E	Fisherbrand* Bags, controlled environment; Bags, polyethylene; 4mil thick; 12L x 10 in. W (30.5 x 25.4cm)	Case of 200	1658
0100292	Fisherbrand* Bags, Sample; Opaque; Sterile; Highly resistant; RNase, DNase, and Pyrogen free; With write-on strip; Safety-tab closures; 3.0 mil thickness; 4.5 x 9 in.; 18 oz. capacity;	Pack of 500	7153
0181810	Fisherbrand* Sample Bag; Heavy-paper; Lined with 1mil polyethylene; Heat sealed; Tin-tie closure; Ideal for collection, shipment or storage of grains; Size: 3-37 x 2.5 x 9.62 in.	Case of 1000	24471
0181811	Fisherbrand* Sample Bag; Heavy-paper; Lined with 1mil polyethylene; Heat sealed; Tin-tie closure; Ideal for collection, shipment or storage of grains; Size: 5 x 3 x 9 in.	Case of 500	15295
0181285A	Fisherbrand* Sample Bag; Disposable; Sterile; Polyethylene; Plain; 3 x 7 in. (8 x 18cm); Capacity, 120mL	Pack of 500	3817
0181285C	Fisherbrand* Sample Bag; Disposable; Sterile; Polyethylene; Plain; 4-1/2 x 9 in. (11 x 23cm); Capacity, 540mL	Pack of 500	5409
14955186	Fisherbrand* Bag; Sample; Polyethylene; Sterile; Clear; With flat-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 18 oz.; Thickness: 3mil; Size: 4-1/2 x 9 in.	Pack of 500	2712
14955187	Fisherbrand* Bag; Sample; Polyethylene; Sterile; Clear; With flat-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 24 oz.; Thickness: 3.5mil; Size: 5-1/2 x 9 in.	Pack of 500	4353
14955188	Fisherbrand* Bag; Sample; Polyethylene; Sterile; Clear; With flat-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 42 oz.; Thickness: 3.5mil; Size: 5-1/2 x 15 in.	Pack of 500	5607
14955189	Fisherbrand* Bag; Sample; Polyethylene; Sterile; Printed; With flat-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 18 oz.; Thickness: 3mil; Size: 4-1/2 x 9 in.	Pack of 500	2996
14955190	Fisherbrand* Bag; Sample; Polyethylene; Sterile; Printed; With flat-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 24 oz.; Thickness: 3.5mil; Size: 5-1/2 x 9 in.	Pack of 500	4797
14955193	Fisherbrand* Bag; Sample; Polyethylene; Sterile; With flat-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 55 oz.; Thickness: 3.5mil; Size: 7 x 12 in.; Clear	Pack of 250	3644
14955194	Fisherbrand* Bag; Sample; Polyethylene; Sterile; With flat-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 55 oz.; Thickness: 3.5mil; Size: 7 x 12 in.; Printed	Pack of 250	4049
14955195	Fisherbrand* Bag; Sample; Polyethylene; Sterile; with flat-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 60 oz.; Thickness: 4mil; Size: 10 x 10 in.; Printed	Pack of 250	5343
14955196	Fisherbrand* Bag; Sample; Polyethylene; Sterile; Clear; With flat-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 100 oz.; Thickness: 4mil; Size: 10 x 14 in.	Pack of 250	5526
0181520	Fisherbrand* Bag, Sample; Printed; Sterile; Polyethylene; Accommodates solid, semi-solid and liquid materials; Round-wire closures; Regular tabs; 2.5mil; 3 x 5 in. (8 x 13cm); 2 oz. capacity	Pack of 500	4622

Part No.	Description	Pack	Price in Rs.
14955175	Fisherbrand* Bag; Sample; Polyethylene; Sterile; Clear; With round-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 2 oz.; Thickness: 2.5mil; Size: 3 x 5 in.	Pack of 500	1195
14955176	Fisherbrand* Bag; Sample; Polyethylene; Sterile; Clear; With round-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 4 oz.; Size: 3 x7 in.; 3.5mil thickness	Pack of 500	1397
14955177	Fisherbrand* Bag; Sample; Polyethylene; Sterile; Clear; With round-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 18 oz.; Thickness: 3mil; Size: 4-1/2 x 9 in.	Pack of 500	2712
14955178	Fisherbrand* Bag; Sample; Polyethylene; Sterile; Clear; With round-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 27 oz.; Thickness: 3mil; Size: 4-1/2 x 12 in.	Pack of 500	3846
14955179	Fisherbrand* Bag; Sample; Polyethylene; Sterile; Clear; With round-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 36 oz.; Thickness: 3mil; Size: 4-1/2 x 15 in.	Pack of 500	4433
14955180	Fisherbrand* Bags, Sample; Polyethylene; Sterile; Clear; With round-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 24 oz.; Size: 5.5 x 9 in.	Pack of 500	4353
14955181	Fisherbrand* Bags, Sample; Polyethylene; Sterile; Clear; With round-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 42 oz.; Thickness: 3.5mil; Size: 5.5 x 15 in.	Pack of 500	5607
14955182	Fisherbrand* Bags, Sample; Polyethylene; Sterile; Printed; With round-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 2 oz.; Thickness: 2.5mil; Size: 3 x 5 in.	Pack of 500	1316
14955183	Fisherbrand* Bags, Sample; Polyethylene; Sterile; Printed; With round-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 4 oz.; Thickness: 4mil; Size: 3 x 7 in.	Pack of 500	1539
14955184	Fisherbrand* Bags, Sample; Polyethylene; Sterile; Printed; With round-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 18 oz.; Thickness: 3mil; Size: 4-1/2 x 9 in.	Pack of 500	2996
14955185	Fisherbrand* Bag; Sample; Polyethylene; Sterile; Printed; With round-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 24 oz.; Thickness: 3.5mil; Size: 5-1/2 x 9 in.	Pack of 500	4797
14955191	Fisherbrand* Bag; Sample; Polyethylene; Sterile; Clear; With round-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 55 oz.; Thickness: 4mil; Size: 7 x 12 in.	Pack of 250	3644
14955192	Fisherbrand* Bag; Sample; Polyethylene; Sterile; Printed; With round-wire closure; Hold solids, semi-solids, liquids; Meet FDA reqts.; Capacity: 55 oz.; Thickness: 4mil; Size: 7 x 12 in.	Pack of 250	3846
0100253	Fisherbrand* Sterile Sampling Bag 100OZ PRT FW	Pack of 250	7054
0100250	Fisherbrand* Sterile Sampling Bag 60OZ CLR FW	Pack of 250	5906
0100290	Fisherbrand* Water Sampling Bag, 2.5 mil thick, 3 x 7 in., 4 oz.,	Case of 100	2973
0100291	Fisherbrand* Water Sampling Bag, 3.0 mil thick, 4.5 x 9 in., 18 oz.,	Case of 100	4817
018161A	Fisherbrand* Bag, Zipper-Seal, Polyethylene, 2 mil thick, 2 x 3 in.	Case of 1000	558
018161B	Fisherbrand* Bag, Zipper-Seal, Polyethylene, 2 mil thick, 3 x 5 in.	Case of 1000	971
018161C	Fisherbrand* Bag, Zipper-Seal, Polyethylene, 2 mil thick, 4 x 6 in.	Case of 1000	867
018161D	Fisherbrand* Bag, Zipper-Seal, Polyethylene, 2 mil thick, 5 x 8 in.	Case of 1000	1385
018161E	Fisherbrand* Bags, Sample; Zipper seal; Transparent polyethylene; Clear, print-free surface; 2mil thick; W x L: 9 x 12 in. (23 x 30cm)	Case of 1000	3243
01816A	Fisherbrand* Bag, Sample, (Zipper) seal, Transparent, Polyethylene, Print-free surface, Not suggested for use with liquids, 4mil thick, 4W x 6 in. Long (10W x 15cmL)	Pack of 50	109
01816B	Fisherbrand* Bag, Sample, (Zipper) seal, Transparent, Polyethylene, Print-free surface, Not suggested for use with liquids, 4mil thick, 6W x 13 in. Long (15W x 33cmL)	Pack of 50	275
01816C	Fisherbrand* Bag, Sample, (Zipper) seal, Transparent, Polyethylene, Print-free surface, Not suggested for use with liquids, 4mil thick, 9W x 18 in. Long (23W x 46cmL)	Pack of 50	630
01816D	Fisherbrand* Bag, Sample; Zipper-seal; Transparent; Polyethylene; Print-free surface; Not suggested for use with liquids; 4mL thick; 5 W x 8 in. Long	Pack of 50	171
01816E	Fisherbrand* Bag, Sample; Zipper seal; Transparent; Polyethylene; Print-free surface; Not suggested for use with liquids; 4mil thick; 9 W x 13 in. L	Pack of 50	514
01816F	Fisherbrand* Bag, Zipper-Seal, Polyethylene, 2 mil thick, 9 x 26 in.	Pack of 50	952
0181314	Fisherbrand* Bags, Sample; Adhesive closure; Tamper evident seal; Non-perforated; Pull away release liner; Saddlepack system; Built-in front and back pouches; 5 x 9 in.	Case of 1000	5463
0181315	Fisherbrand* Bags; Adhesive Closure; Tamper Evident; Non-perforated; Pull away release liner; Built-in front and back pouches; 6 x 10 in.	Case of 1000	6402