

Part No.	Description	Pack	Price in Rs.
035415	Fisherbrand* Brush, Heavy Duty Bristle Beaker Large Brush, 16 in.	Each	274
FBA03562	Fisherbrand* Brush, Cone; Imhoff Settling; Double bristle fan tip flares for increased scrubbing area; Tapered brush head; Smooth wooden handle; Natural and nylon bristles Overall Length: 26 in. (66cm); Bristle L x Dia.: 91/2 in. (24cm)	Each	344
FBA03661	Fisherbrand* Brush, Balances; Slightly pointed tip and full body tapering toward handle; Bristle: 1-1/2 x 3/4 in.; Length: 9 in.	Each	384
FBA03555	Fisherbrand* Brush, Beaker; End tuft for cleaning corners; Wooden handle; Two rows of bristles provide slightly oval shape; Bristles: 6 x 3 in.; Length: 16 in.	Each	208
03621A	Fisherbrand* Brush, Nessler Tubes; Low-Form; double-fan tip; Full bristles; galvanized wire handle; Overall Length: 12 in.; Bristles: 4 L x 1.75 in. Dia.; Small and Large Cylinder;	Each	96
03621B	Fisherbrand* Brush, Nessler Tubes; Radial, double-fan tip; For low-form tubes; 20 in. L Overall; Bristle dia. 2.38 in. x 6 in. L	Each	208
03621C	Fisherbrand* Brush, Nessler Tubes; Radial, double-fan tip; For low-form tubes; 30 in. L Overall; Bristles: 6 L x 2.38 in. Dia.	Pack of 12	3072
FBA03622	Fisherbrand* Brush, Nessler Tubes; Radial, double-fan tip; For high-form tubes; 19 in. L; Bristles: 4 L x 1 in. Dia.	Each	104
FBA03590	Fisherbrand* Brushes, Test Tube; Nylon Tufted Brush; Black, nylon tuft, 1 in. (2.5cm); Bleached white hog bristles; Galvanized wire handle; Overall Length: 10 in. (25cm); Bristle L x Dia. : 3 x 13/8 in. (7.6 x 3.5cm)	Pack of 12	576
FBA03670	Fisherbrand* Brush, Round Camel's Hair; Soft brown hair set into metal ferrule mounted on thin wooden handle; For dusting or painting; Overall Length: 7-3/4 in. (19cm); Bristle L x Dia. :13/16 x 3/16 in. (2.1 x 0.5cm)	Pack of 12	480
FB70066	Fisherbrand* Sieve brush double ended brass 30mm x 135mm	Each	552
FB70065	Fisherbrand* Sieve brush double ended nylon and bristle 20mm x 230mm	Each	610
FBA03636	Fisherbrand* Brush, Capillary tube; Soft Nylon Bristle; For cleaning funnel stems and capillary tubes; Bristle: 2-1/2 x 1/2 in.; Length: 12 in.	Pack of 12	403
035725	Fisherbrand* Brush, Test tube; Stiff, white nylon bristles; Radial tufted end; Stainless-steel handle is rustproof, corrosion resistant; Overall Length: 11 in.; Bristle L x D: 4 x 1-3/8 in.	Pack of 12	1728
FBA03578	Fisherbrand* Brush, Test tube; Stiff, white nylon bristles; Radial tufted end; Stainless-steel handle is rustproof, corrosion resistant; Bristles: 3/4 x 3-1/2 in.; Length: 9 in.	Pack of 12	960
034051	Fisherbrand* StaticAway* Brush, Anti-Static; Dissipates Static and Whisk Harmful Dust, Lint, Dirt	Each	2272
FBA03635	Fisherbrand* Brushes; Stiff Bristled; For Clean tubes to 7/16 in. dia. (1.1cm); Overall Length: 12 in. (30cm); Bristles: 3L x 1/2 in. Dia. (7.6 x 1.3cm)	Pack of 12	480
FBA03682	Fisherbrand* Stiff Palmyra Brush with Handle Overall Length: 8 3/8 in.	Each	613
FBA03680	Fisherbrand* Brush, Table dusting; Soft-bristle brush; Handle, 5 in.L (13cm); For dusting tables and lab benches; Overall Length: 125/8 in. (32cm); Bristle, L x W : 71/2 x 2 in. (19 x 5cm)	Each	807
FBA03625	Fisherbrand* Brush Transfer Pipet Stiff Nylon Bristles 17 in.	Each	61
FBA03536	Fisherbrand* Brush, Beaker; White; Nylon; Bristle; Beaker Brush; Wooden handle; Chemical-resistant; Overall Length 16 in. (41cm); Bristle L x Dia. 5 x 27/8 in. (13 x 7cm)	Each	274
60180778	Thermo Scientific* Injector Cleaning Brush	Each	596

Buckets

Part No.	Description	Pack	Price in Rs.
FBA03438	Fisherbrand* Pail; Heavy-duty HDPE; Large spout, splash shield, and handle; Graduated in 1 gal. and 1L increments; Capacity, 4 gal.	Each	9219