

## Burets

Part No.	Description	Pack	Price in Rs.
<b>Glass General Purpose Burets</b>			
03700A	Fisherbrand* Buret, Manual; With blue markings; With straight standard taper glass stopcock; Capacity: 10mL	Each	1981
03700B	Fisherbrand* Buret; Manual; With blue markings; Capacity: 25mL	Each	1498
03700C	Fisherbrand* Buret; Manual, With blue markings, Capacity: 50mL	Each	2378
03700D	Fisherbrand* Buret; Manual, With blue markings, Capacity: 100mL	Each	2496
03740A	Fisherbrand* Buret, Schellbach; Straight standard taper glass stopcock; Accentuate meniscus; Accuracy established w/inthe limits of error allowed by NIST; Capacity: 25mL	Each	3738
03740B	Fisherbrand* Buret, Schellbach; Straight standard taper glass stopcock; Accentuate meniscus; Accuracy established w/inthe limits of error allowed by NIST; Capacity: 50mL	Each	3330
03740C	Fisherbrand* Buret, Schellbach; Straight standard taper glass stopcock; Accentuate meniscus; Accuracy established w/inthe limits of error allowed by NIST; Capacity: 100mL	Each	4179
03750A	Fisherbrand* Buret, Schellbach; W/3-way stopcock; Automatic zeroing; Black graduations; Calibrated to a tip at top; Calibration accuracy to ASTM E 287; Capacity: 25mL	Each	11218
03750B	Fisherbrand* Buret, Schellbach; W/3-way stopcock; Automatic zeroing; Black graduations; Calibrated to a tip at top; Calibration accuracy to ASTM E 287; Capacity: 50mL	Each	11749
03750C	Fisherbrand* Buret, Schellbach; W/3-way stopcock; Automatic zeroing; Black graduations; Calibrated to a tip at top; Calibration accuracy to ASTM E 287; Capacity: 100mL	Each	13779
<b>Glass Automatic Burets</b>			
038421A	Fisherbrand* Buret; With clear glass reservoir; Reservoir Capacity: 500mL; Capacity: 10mL	Each	4839
038421B	Fisherbrand* Buret, Automatic; With clear glass reservoir; Reservoir Capacity: 500mL; Capacity: 25mL; Subdivision: 0.1mL	Each	4852
038422B	Fisherbrand* Buret, Automatic; With amber glass reservoir; Reservoir Capacity: 500mL; Capacity: 25mL; Subdivision: 0.1mL	Each	4951
038423A	Fisherbrand* Buret, Automatic; with Std. Taper Stopcocks and Storage Bottles Replacement Burets; Capacity: 10mL	Each	2938
038423B	Fisherbrand* Automatic Buret with Std. Taper Stopcocks and Storage Bottles Replacement Burets Capacity: 25mL	Each	3029
038431A	Fisherbrand* Buret, Automatic; Includes: Glass buret w/central filling tube, 500mL clear glass reservoir and delivery tip, rubber stopper and bulb, glass air tube and bead; Pinch valve; Cap.: 10mL; Sub.: 0.10mL	Each	4012
038431B	Fisherbrand* Buret, Automatic; Includes: Glass buret w/central filling tube, 500mL clear glass reservoir and delivery tip, rubber stopper and bulb, glass air tube and bead; Pinch valve; Cap.: 25mL; Sub.: 0.10mL	Each	4839
038432A	Fisherbrand* Buret, Automatic; Includes: Glass buret w/central filling tube, 500mL amber glass reservoir and delivery tip, rubber stopper and bulb, glass air tube and bead; Pinch valve; Cap.: 10mL; Sub.: 0.10mL	Each	4130
038432B	Fisherbrand* Buret, Automatic; Includes: Glass buret w/central filling tube, 500mL amber glass reservoir and delivery tip, rubber stopper and bulb, glass air tube and bead; Pinch valve; Cap.: 25mL; Sub.: 0.10mL	Each	4396
038433A	Fisherbrand* Buret, Replacement; For 03-843-1 and 03-843-2 series; Clear; Capacity: 10mL	Each	3413
038433B	Fisherbrand* Buret, Replacement; For 03-843-1 and 03-843-2 series; Clear; Capacity: 25mL	Each	2971
03845AAT	Fisherbrand* Buret, Automatic; Class A; For petroleum neutralization; With stopcock, reagent bottle, and aspirator bulb; Reservoir capacity: 1000mL; Capacity: 5mL	Each	22632
03845AT	Fisherbrand* Buret, Automatic; Class A; For petroleum neutralization; With stopcock, reagent bottle, and aspirator bulb; Reservoir capacity: 1000mL; Capacity: 10mL	Each	22568
03845BBT	Fisherbrand* Buret, Automatic; Class A; For petroleum neutralization; With stopcock, reagent bottle, and aspirator bulb; Reservoir capacity: 1000mL; Capacity: 25mL	Each	25002
03845BT	Fisherbrand* Buret, Automatic; Class A; For petroleum neutralization; With stopcock, reagent bottle, and aspirator bulb; Reservoir capacity: 1000mL; Capacity: 50mL	Each	27303
03846B	Fisherbrand* Auto-Zero Class A Buret; Acid clear; Capacity, 50mL; 2mm glass stopcock; With reservoir and white graduations	Each	20659
038471A	Fisherbrand* Machlett Automatic Burets Reservoir Capacity: 500mL; Capacity: 2mL	Each	18830
038471B	Fisherbrand* Buret, Automatic Machlett; PTFE PTFE-plug stopcock; Subdivision 0.02mL; Reservoir Capacity: 1000mL; Buret Capacity: 5mL	Each	19675
038471C	Fisherbrand* Buret, Matchlett Automatic; Class A; Capacity: 10mL; Reservoir Capacity: 1000mL; Subdivision: 0.05mL; Tolerance; +/- 0.02mL	Each	21317
038471D	Fisherbrand* Buret, Matchlett Automatic; Class A; Capacity: 25mL; Reservoir Capacity: 2000mL; Subdivision: 0.05mL; Tolerance; +/- 0.02mL	Each	22995
038471E	Fisherbrand* Buret, Automatic Machlett; PTFE-plug stopcock; Subdivision 0.1mL; Reservoir Capacity: 2000mL; Buret Capacity: 50mL	Each	24819
<b>Dispensing Burets</b>			
03790A	Fisherbrand* Buret, Dispensing; With straight glass stopcock; Capacity: 250mL; Subdivisions: 1mL	Each	9336
03846AT	Fisherbrand* Buret; Acid-Dispensing; With PTFE stopcock, 1000mL reservoir and aspirator bulb; Capacity: 10mL	Each	20113

Part No.	Description	Pack	Price in Rs.
03846BBT	Fisherbrand* Buret; Acid-Dispensing; With PTFE stopcock, 1000mL reservoir and aspirator bulb; Capacity: 25mL	Each	21350
03846BT	Fisherbrand* Buret; Acid-Dispensing; With PTFE stopcock, 1000mL reservoir and aspirator bulb; Capacity: 50mL	Each	22933
<b>Digital Burets</b>			
92003	Fisher Scientific* Qualigens* DIGITRATE, Titration Burette 30.0 ml, Increments - 0.01 ml	Each	45000
92004	Fisher Scientific* Qualigens* DIGITRATE, Titration Burette 50.0 ml, Increments - 0.01 ml	Each	45000
<b>Buret Accessories</b>			
FBA03872	Fisherbrand* Meniscus Reader, With Vinylite reader and transparent front strap with snap to hold reader in place	Each	472

## Burners

Part No.	Description	Pack	Price in Rs.
<b>Gas Burners</b>			
03902Q	Fisher Scientific* accuFlame* Bunsen Burner; Stainless steel grid, chrome-plated brass body, Zn-Al base; Natural gas; Heat output 2930w/10,000BTU/hr.; Grid dia. 3.8cm; Ht. 22cm; No-tip base; For use with 1/4 in. I.D. gas tubing, not inc.	Each	2506
03917Q	Fisher Scientific* Economy Burners Heat Value: 800-1200BTU/ft.3 (29.8-44.7MJ/m3); Type: Natural Gas Type	Each	1101
04067AQ	Fisher Scientific* Small-Size Burners Heat Value: 150-600 BTU/ft.3(5.6-22.4MJ/m3); Type: Artificial Gas Type	Each	1506
04067NQ	Fisher Scientific* Small-Size Burners Heat Value: 800-1200BTU/ft.3 (29.8-44.7MJ/m3); Type: Natural Gas Type	Each	1247
03962PQ	Fisher Scientific* Adjustable-Flame Burners Heat Value: 2400-3300 BTU/ft.3 (89.4-122.9MJ/m3); Type: LP Gas Type	Each	1798
03962Q	Fisher Scientific* Adjustable-Flame Bunsen Burner; Pittsburgh Universal straight tube; Heat output 5000 BTU/hr.; 13 x 165mm; 800-1200BTU/cu. ft. (29.8-44.7MJ/m3); Natural Gas	Each	1815
03900Q	Fisher Scientific* Standard Grid-Top Adjustable-Flame Burners Heat Value: 150-600 BTU/ft.3 (5.6-22.4MJ/m3); Type: Artificial Gas Type	Each	3156
03902PQ	Fisher Scientific* Standard Grid-Top Adjustable-Flame Burners Burners, Standard Grid-Top; Fisher Scientific; Adjustable-Flame; Channel flame through 100 separate openings; Rustproof chromium-plated brass body; Grid Outside Dia.: 1-1/2 in. (3.8cm); Heat Value 2400 - 3300 BTU/ft3 ((89.4-122.9MJ/m3); Height: 8-5/8 in.	Each	3005
<b>Burner Accessories</b>			
03908Q	Fisher* accuFlame* Burner Cap and Grid for Natural Gas Burner 1-1/2 in	Each	799
00146Q	Fisher* Burner Halo Support for Burners 2-1/2 in	Each	909
S41893D	Fisherbrand* Matches, Safety; Strike-on-box; Approx. 28 matches per box	Pk of 10	116

## Calculators

Part No.	Description	Pack	Price in Rs.
S40678	Fisher Scientific* Conversion calculator, Provides over 100 conversions, Plus a full-function 8 digit calculator	Each	3600
1209923	Fisherbrand* Calculator, 12-digit; 1/2 in. (1.3cm) display; Features memory, percent, floating decimal point, ratios, and grand totals; Includes solar panels, silver-oxide battery, and ABS plastic case; 4 1/2L x 2 1/2W x 5/8 in.H (11.4 x 6.4 x 1.6cm); 3 o	Each	4004
14648100	Fisherbrand* Calculator; Clipboard; With Retractable Measuring Tape; Large display; Solar-powered calculator; Clipboard: 14 x 9 in.; Tape measures to 60 in. (150cm)	Each	3909
14648230	Fisherbrand* Calculator; Clipboard; With Retractable Measuring Tape; Large display; Solar-powered calculator; Clipboard: 14 x 9 in.; Tape measures to 60 in. (150cm)	Each	3909
1464850	Fisherbrand* Calculator, Conversion; 126 types of metric conversions; Converts length, area, volume, weight, liquid, energy, and temperature; Performs English (Imperial) conversions, and has all routine calculator functions; 1/2 in. digits; 8-digit LCD;	Each	2656
FBA12099	Fisherbrand* Calculator, Metric-Conversion; 8 digit 1/2-inch-high LCD; Converts length, area, volume, weight, liquid, energy, and temperature; Performs English (Imperial) conversions and all routine calculator functions; 3 3/4L x 3W 1/2 in.H (9.5 x 7.6 x	Each	5810
1209921	Fisherbrand* Calculator, Pocket-sized, Conversion, 4-1/2 L x 2-5/8 W x 5/16 in H	Each	1639
1209924	Fisherbrand* Calculator, Solar; Desktop; 8-digit, 3/4 in. display; ABS plastic case; Protective feet; 6 1/2L x 4 1/2W x 3/16 in.H; 4 oz.	Each	1458
1209925	Fisherbrand* Calculator, Solar; Desktop; 12-digit; Calculates up to 999,999,999,999 and down to .000000000001; 3/4 in. digits, ABS plastic case; Protective feet; 6 1/2L x 4 1/2W x 3/16 in.H; 4 oz.	Each	1895
1209920	Fisherbrand* Calculator, Big digit, Solar-powered, 4-1/4 L x 2-3/8 W x 1/4 in H	Each	954
1209912	Fisherbrand* Universal Converter/Timer Weight: 7/8 oz.; Dimensions: 2 1/4L x 3 3/8W x 1/4 in.H	Each	2256