

Counters

Part No.	Description	Pack	Price in Rs.
Frequency Counters			
1464977	Fisherbrand* Pen, counter; Combination marker/digital counter; Touch, sound, and count occur simultaneously; Counter reads to 99,999; Fail-safe reset	Each	1435
1464978	Fisherbrand* Counter, Digital; Key-Chain/Wrist-Strap/Belt-Clip; LCD; Recessed failsafe reset; Water-resistant; Shock-resistant; O-ring-sealed ABS plastic case; Button-click for count confirmation	Each	983
Parts Counters			
FB57942	Fisherbrand* Counter with total window each key records 999 250mm x 76mm x 56mm mechanical 6 key	Each	22032
FB57944	Fisherbrand* Counter with total window each key records 999 320mm x 76mm x 56mm mechanical 9 key	Each	25245
079052	Fisherbrand* Counter, Digital Reads to 99,999; Water- and shock-resistant; Easy-to-read 1/4 in. LCD; Solid-state construction; Battery power for two years' continuous service	Each	1480
079056	Fisherbrand* Counter, Hand Tally; Tallies to 9999; Resets to 0; Pocket size; Chrome finish	Each	539

Crucibles

Part No.	Description	Pack	Price in Rs.
Glass Crucibles			
104904	Fisherbrand* Crucibles, Glass; Quartz, Self-sealing; O.D. x H: 7/8 x 1-5/16 in. (22.5 x 34mm)	Case of 4	28133
Metal Crucibles			
13812120	Fisherbrand* Crucible, Nickel Chromium; One piece; For avolatile matter tests; Withstands temps to 1000degC; 13mL capacity; 26mm ID x 34mm H; Without cover	Each	2800
13812121	Fisherbrand* Crucible Cover; For Nickel Chromium Crucible; 26mm ID x 11mm H	Each	1040
13812126	Fisherbrand* Crucible, Nickel; One piece; For diluting alkalis; Withstands temps to 800degC; 20mL capacity; 38mm ID x 32mm H; Without cover	Each	1360
13812127	Fisherbrand* Crucible, Nickel; One piece; For diluting alkalis; Withstands temps to 800degC; 30mL capacity; 44mm ID x 44mm H; Without cover	Each	1600
13812128	Fisherbrand* Crucible, Nickel; One piece; For diluting alkalis; Withstands temps to 800degC; 50mL capacity; 44mm ID x 51mm H; Without cover	Each	2080
13812129	Fisherbrand* Crucible, Nickel; One piece; For diluting alkalis; Withstands temps to 800degC; 75mL capacity; 57mm ID x 57mm H; Without cover	Each	2560
13812130	Fisherbrand* Crucible, Nickel; One piece; For diluting alkalis; Withstands temps to 800degC; 100mL capacity; 64mm ID x 57mm H; Without cover	Each	3200
13812131	Fisherbrand* Crucible, Nickel; One piece; For diluting alkalis; Withstands temps to 800degC; 250mL capacity; 84mm ID x 78mm H; Without cover	Each	5761
13812132	Fisherbrand* Crucible Cover; Nickel; For 20mL Nickel Crucible	Each	640
13812133	Fisherbrand* Crucible Cover; Nickel; For 30mL Nickel Crucible	Each	760
13812134	Fisherbrand* Crucible Cover; Nickel; For 50mL Nickel Crucible	Each	880
13812135	Fisherbrand* Crucible Cover; Nickel; For 75mL Nickel Crucible	Each	1240
13812136	Fisherbrand* Crucible Cover; Nickel; For 100mL Nickel Crucible	Each	1424
13812150	Fisherbrand* Crucibles, Metal; Cover; For 20mL steel crucible	Case of 10	5601
13812151	Fisherbrand* Crucible, Steel; Low-form one-piece sheet steel; For alkalis and peroxide fusions; Withstands temps to 500degC; 20mL capacity; 36mm ID x 26mm H; Without cover;	Pack of 10	5601
13812152	Fisherbrand* Crucible, Steel; Low-form one-piece sheet steel; For alkalis and peroxide fusions; Withstands temps to 500degC; 100mL capacity; 56mm ID x 48mm H; Without cover;	Each	1040
13812153	Fisherbrand* Crucible, Steel; Low-form one-piece sheet steel; For alkalis and peroxide fusions; Withstands temps to 500degC; 200mL capacity; 68mm ID x 60mm H; Without cover;	Each	1360

Part No.	Description	Pack	Price in Rs.
13812154	Fisherbrand* Crucible, Steel; Low-form one-piece sheet steel; For alkalis and peroxide fusions; Withstands temps to 500degC; 400mL capacity; 94 ID x 65mm H; Without cover;	Pack of 10	17602
13812122	Fisherbrand* Crucible, Zirconium; One piece polished; Withstands up to 100 fusions (2-5 min. between 460degC and 900degC); 25mL capacity; 44mm ID x 22mm H)	Each	5601
13812123	Fisherbrand* Crucible, Zirconium; One piece polished; Withstands up to 100 fusions (2-5 min. between 460degC and 900degC); 35mL capacity; 44mm ID x 29mm H)	Each	5761
13812124	Fisherbrand* Crucible, Zirconium; One piece polished; Withstands up to 100 fusions (2-5 min. between 460degC and 900degC); 45mL capacity; 44mm ID x 34mm H)	Each	5921
13812125	Fisherbrand* Crucible, Zirconium; One piece polished; Withstands up to 100 fusions (2-5 min. between 460degC and 900degC); 55mL capacity; 44mm ID x 41mm H)	Each	6161
PTFE Crucibles			
02617131	Fisherbrand* Crucibles, PTFE; 5mL; Diameter x Height: 1-1/16 x 5/8 in.; Withstands temperatures to 280 deg. C (536 deg. F)	Each	826
02617132	Fisherbrand* Crucibles, PTFE; 25mL; Diameter x Height: 1-3/4 x 1 in.; Withstands temperatures to 280 deg. C (536 deg. F)	Each	1202
02617133	Fisherbrand* Crucibles, PTFE; 75mL; Diameter x Height: 2-1/4 x 1-5/8 in.; Withstands temperatures to 280 deg. C (536 deg. F)	Each	1975
02617134	Fisherbrand* Crucibles, PTFE; 100mL; Diameter x Height: 2-13/16 x 1-7/8 in.; Withstands temperatures to 280 deg. C (536 deg. F)	Each	2290

Cryogenics

Part No.	Description	Pack	Price in Rs.
Cryogenic Storage			
FB57938	Fisherbrand* 10 x 10 array for 1.2mL/2mL vials, includes picking tool	Pack of 4	4113
FB57935	Fisherbrand* 5 x 5 array for 1.2mL/2mL vials	Pack of 8	5095
FB57936	Fisherbrand* 9 x 9 array for 1.2mL/2mL vials, includes picking tool	Pack of 4	4163
FB57937	Fisherbrand* 9 x 9 array tall boxes for 3mL/5mL vials	Pack of 5	6018
03391200	Fisherbrand* Cryogenic Storage Box; Polypropylene; 100 places; For 0.5/1.5/2.0mL tubes/cryo vials up to 5.35cm high; Temp. range -80degree to +121degreeC; Markable; Dual alphanumeric grid on lid; Stackable; 14L x 14W x 6cm H; Natural	Pack of 5	3000
03391201	Fisherbrand* Cryogenic Storage Box; Polypropylene; 100 places; For 0.5/1.5/2.0mL tubes/cryo vials up to 5.35cm high; Temp. range -80degree to +121degreeC; Markable; Dual alphanumeric grid on lid; Stackable; 14L x 14W x 6cm H; Blue	Pack of 5	3000
03391202	Fisherbrand* Cryogenic Storage Box; Polypropylene; 100 places; For 0.5/1.5/2.0mL tubes/cryo vials up to 5.35cm high; Temp. range -80degree to +121degreeC; Markable; Dual alphanumeric grid on lid; Stackable; 14L x 14W x 6cm H; Green	Pack of 5	3000
03391203	Fisherbrand* Cryogenic Storage Box; Polypropylene; 100 places; For 0.5/1.5/2.0mL tubes/cryo vials up to 5.35cm high; Temp. range -80degree to +121degreeC; Markable; Dual alphanumeric grid on lid; Stackable; 14L x 14W x 6cm H; Purple	Pack of 5	3000
03391204	Fisherbrand* Cryogenic Storage Box; Polypropylene; 100 places; For 0.5/1.5/2.0mL tubes/cryo vials up to 5.35cm high; Temp. range -80degree to +121degreeC; Markable; Dual alphanumeric grid on lid; Stackable; 14L x 14W x 6cm H; Orange	Pack of 5	3000
03391205	Fisherbrand* Cryogenic Storage Box; Polypropylene; 100 places; For 0.5/1.5/2.0mL tubes/cryo vials up to 5.35cm high; Temp. range -80degree to +121degreeC; Markable; Dual alphanumeric grid on lid; Stackable; 14L x 14W x 6cm H; Assorted	Pack of 5	3000
03395400	Fisherbrand* Cryogenic Storage Box; Polypropylene; 81 places; For 0.5/1.5/2.0mL tubes/cryo vials up to 5.35cm high; Temp. range -80degree to +121degreeC; Markable; Dual alphanumeric grid on lid; Stackable; 12.7L x 12.7W x 5.1cm H; Natural	Pack of 5	2352
03395401	Fisherbrand* Cryogenic Storage Box; Polypropylene; 81 places; For 0.5/1.5/2.0mL tubes/cryo vials up to 5.35cm high; Temp. range -80degree to +121degreeC; Markable; Dual alphanumeric grid on lid; Stackable; 12.7L x 12.7W x 5.1cm H; Assorted	Pack of 5	2352
03395402	Fisherbrand* Cryogenic Storage Box; Polypropylene; 81 places; For 0.5/1.5/2.0mL tubes/cryo vials up to 5.35cm high; Temp. range -80degree to +121degreeC; Markable; Dual alphanumeric grid on lid; Stackable; 12.7L x 12.7W x 5.1cm H; Blue	Pack of 5	2352
03395403	Fisherbrand* Cryogenic Storage Box; Polypropylene; 81 places; For 0.5/1.5/2.0mL tubes/cryo vials up to 5.35cm high; Temp. range -80degree to +121degreeC; Markable; Dual alphanumeric grid on lid; Stackable; 12.7L x 12.7W x 5.1cm H; Purple	Pack of 5	2352
03395404	Fisherbrand* Cryogenic Storage Box; Polypropylene; 81 places; For 0.5/1.5/2.0mL tubes/cryo vials up to 5.35cm high; Temp. range -80degree to +121degreeC; Markable; Dual alphanumeric grid on lid; Stackable; 12.7L x 12.7W x 5.1cm H; Orange	Pack of 5	2352
14826340	Fisherbrand* Cryogenic Storage Box; Polypropylene; 81 places; For 0.5/1.5/2.0mL tubes/cryo vials up to 5.35cm high; Temp. range -80degree to +121degreeC; Markable; Dual alphanumeric grid on lid; Stackable; 12.7L x 12.7W x 5.1cm H; Green	Pack of 5	2352
1050028	Fisherbrand* Rack; Cryogenic vial storage box; White Polycarbonate; Cover printed with 9 x 9 grid; Autoclavable; Stackable; For internally and externally threaded vials; For 1.2 to 2.0mL vials; Dimensions: 5-1/4L x 5-1/4W x 2 in.H	Pack of 4	2601