

## Cylinders

Part No.	Description	Pack	Price in Rs.
<b>Glass Graduated Cylinders for Mixing</b>			
08565A	Fisherbrand* Cylinders, Mixing; Calibrated TC; Dual White Scale; 10mL; Subdiv: 0.2mL; w/stopper	Each	1949
08565B	Fisherbrand* Cylinders, Mixing; Calibrated TC; Dual White Scale; 25mL; Subdiv: 0.5mL; w/stopper	Each	2162
08565C	Fisherbrand* Cylinders, Mixing; Calibrated TC; Dual White Scale; 50mL; Subdiv: 1mL; w/stopper	Each	2768
08565D	Fisherbrand* Cylinders, Mixing; Calibrated TC; Dual White Scale; 100mL; Subdiv: 1mL; w/stopper	Each	3203
08565E	Fisherbrand* Cylinders, Mixing; Calibrated TC; Dual White Scale; 250mL; Subdiv: 2mL; w/stopper	Each	5385
08565F	Fisherbrand* Cylinders, Mixing; Calibrated TC; Dual White Scale; 500mL; Subdiv: 5mL; w/stopper	Each	6271
<b>Glass Graduated Cylinders Rated To Contain</b>			
08559AC	Fisherbrand* Cylinders, Graduated; Class A; w/traceability/calib.; TC; White scale; 10mL; Subdiv: 0.1mL	Each	9614
08559BC	Fisherbrand* Cylinders, Graduated; Class A; w/traceability/calib.; TC; White scale; 25mL; Subdiv: 0.2mL	Each	9516
08559CC	Fisherbrand* Cylinders, Graduated; Class A; w/traceability/calib.; TC; White scale; 50mL; Subdiv: 1mL	Each	9457
08559DC	Fisherbrand* Cylinders, Graduated; Class A; w/traceability/calib.; TC; White scale; 100mL; Subdiv: 1mL	Case of 6	59430
08559EC	Fisherbrand* Cylinders, Graduated; Class A; w/traceability/calib.; TC; White scale; 250mL; Subdiv: 2mL	Each	12272
08559FC	Fisherbrand* Cylinders, Graduated; Class A; w/traceability/calib.; TC; White scale; 500mL; Subdiv: 5mL	Each	13205
08559GC	Fisherbrand* Cylinders, Graduated; Class A; w/traceability/calib.; TC; White scale; 1000mL; Subdiv: 10mL	Case of 2	29708
08559HC	Fisherbrand* Cylinders, Graduated; Class A; w/traceability/calib.; TC; White scale; 2000mL; Subdiv: 20mL	Each	20145
08559A	Fisherbrand* Cylinder, Graduated; Serialized Class A; w/traceability; TC; Borosilicate glass; White scale; Capacity: 10mL; Subdiv: 0.1mL; Tolerance: +/-0.10mL; Approx. height: 140mm	Each	2712
08559B	Fisherbrand* Cylinder, Graduated; Serialized Class A; w/traceability; TC; Borosilicate glass; White scale; Capacity: 25mL; Subdiv: 0.5mL; Tolerance: +/-0.17mL; Approx. height: 140mm	Each	2514
08559C	Fisherbrand* Cylinder, Graduated; Serialized Class A; w/traceability; TC; Borosilicate glass; White scale; Capacity: 50mL; Subdiv: 1mL; Tolerance: +/-0.25mL; Approx. height: 190mm	Each	2674
08559D	Fisherbrand* Cylinder, Graduated; Serialized Class A Glass; Calibrated To Contain; Supplied with Certificate of Traceability; White scale graduations; Capacity: 100mL; Subdivisions: 1mL	Each	2766
08559E	Fisherbrand* Cylinder; Serialized class A; With Certificate of Traceability; Borosilicate glass; Graduated; Calibrated to contain; Capacity: 250mL; Subdivisions: 2mL; Tolerance: +/-1mL; Approx. height: 330mm	Each	4430
08559F	Fisherbrand* Cylinder; Serialized class A; With Certificate of Traceability; Borosilicate glass; Graduated; Calibrated to contain; Capacity: 500mL; Subdivisions: 5mL; Tolerance: +/-2mL; Approx. height: 375mm	Each	5235
08559G	Fisherbrand* Cylinders, Graduated; Class A; Serialized; With certificate of traceability; Calibrated To Contain; White scale; Molded hexagonal base; Capacity: 1000mL; Subdiv. 10mL; Tolerance: +/-2.0mL; Approx. Height: 460mm	Each	6725
08559H	Fisherbrand* Cylinder; Serialized class A; With Certificate of Traceability; Borosilicate glass; Graduated; Calibrated to contain; Capacity: 2000mL; Subdivisions: 20mL; Tolerance: +/-6mL; Approx. height: 520mm	Each	13020
085571A	Fisherbrand* Cylinder, Graduated; Class B Glass; Calibrated To Contain; White-ceramic graduations; Handy pouring spouts; Capacity: 10mL; Subdivisions: 0.2mL	Pk of 2	2012
085571B	Fisherbrand* Cylinder, Graduated; Class B Glass; Calibrated To Contain; White-ceramic graduations; Handy pouring spouts; Capacity: 25mL; Subdivisions: 0.5mL	Each	1328
085571C	Fisherbrand* Cylinder, Graduated; Class B Glass; Calibrated To Contain; White-ceramic graduations; Handy pouring spouts; Capacity: 50mL; Subdivisions: 1mL	Each	1384
085571D	Fisherbrand* Cylinder, Graduated; Class B Glass; Calibrated To Contain; White-ceramic graduations; Handy pouring spouts; Capacity: 100mL; Subdivisions: 1mL	Each	1612
085571E	Fisherbrand* Cylinder, Graduated; Class B Glass; Calibrated To Contain; White-ceramic graduations; Handy pouring spouts; Capacity: 250mL; Subdivisions: 2mL	Each	2718

Part No.	Description	Pack	Price in Rs.
085571F	Fisherbrand* Cylinder, Graduated; Class B Glass; Calibrated To Contain; White-ceramic graduations; Handy pouring spouts; Capacity: 500mL; Subdivisions: 5mL	Each	3479
085571G	Fisherbrand* Cylinder, Graduated; Class B Glass; Calibrated To Contain; White-ceramic graduations; Handy pouring spouts; Capacity: 1000mL; Subdivisions: 10mL	Each	5024
085571H	Fisherbrand* Cylinder, Graduated; Class B Glass; Calibrated To Contain; White-ceramic graduations; Handy pouring spouts; Capacity: 2000mL; Subdivisions: 20mL	Each	8589
08563B	Fisherbrand* Cylinders, Graduated; Class B Glass; TC; White English Scale; 4 oz.; Subdiv: 1/16 oz.	Each	1813
08563C	Fisherbrand* Cylinders, Graduated; Class B Glass; TC; White English Scale; 8 oz.; Subdiv: 1/8 oz.	Each	2723
08563D	Fisherbrand* Cylinders, Graduated; Class B Glass; TC; White English Scale; 16 oz.; Subdiv: 1/4 oz.	Each	3176
08563E	Fisherbrand* Cylinders, Graduated; Class B Glass; TC; White English Scale; 32 oz.; Subdiv: 1/2 oz.	Each	4764
<b>Glass Graduated Cylinders Rated To Deliver</b>			
08557AC	Fisherbrand* Cylinders, Graduated; Class A; w/traceability/calib.; TD; White scale; 10mL; Subdiv: 0.1mL	Each	9223
08557BC	Fisherbrand* Cylinders, Graduated; Class A; w/traceability/calib.; TD; White scale; 25mL; Subdiv: 0.2mL	Each	9049
08557CC	Fisherbrand* Cylinders, Graduated; Class A; w/traceability/calib.; TD; White scale; 50mL; Subdiv: 1mL	Each	9070
08557DC	Fisherbrand* Cylinders, Graduated; Class A; w/traceability/calib.; TD; White scale; 100mL; Subdiv: 1mL	Each	9476
08557EC	Fisherbrand* Cylinders, Graduated; Class A; w/traceability/calib.; TD; White scale; 250mL; Subdiv: 2mL	Each	11748
08557FC	Fisherbrand* Cylinders, Graduated; Class A; w/traceability/calib.; TD; White scale; 500mL; Subdiv: 5mL	Each	12623
08557HC	Fisherbrand* Cylinders, Graduated; Class A; w/traceability/calib.; TD; White scale; 2000mL; Subdiv: 20mL	Each	20145
0973031	Fisherbrand* Cylinders, Graduated; Class A; w/traceability/cert.; TD; White scale; 1000mL; Subdiv: 10mL	Each	14488
08557A	Fisherbrand* Cylinders, Graduated; Class A; Serialized; With certificate of traceability; Calibrated To Deliver; White scale; Molded hexagonal base; Capacity: 10mL; Subdiv. 0.1mL; Tolerance: +/- .10mL; Approx. Height: 175mm	Each	2671
08557B	Fisherbrand* Cylinders, Class A; Glass; Serialized; Calibrated To Deliver; With Certificate of Traceability; White scale; Molded hexagonal base; Capac.: 25mL; Subdiv.: 0.2mL; Tol.: +/-0.17mL; Approx. Height: 200mm	Each	2439
08557C	Fisherbrand* Cylinders, Class A; Glass; Serialized; Calibrated To Deliver; With Certificate of Traceability; White scale; Molded hexagonal base; Capac.: 50mL; Subdiv.: 1mL; Tol.: +/-0.25mL; Approx. Height: 225mm	Each	2528
08557D	Fisherbrand* Cylinder, Graduated; Serialized Class A Glass; Calibrated To Deliver; With Certificate of Traceability; White scale graduations; Capacity: 100mL; Subdivisions: 1mL	Each	2943
08557E	Fisherbrand* Cylinder, Graduated; Serialized Class A Glass; Calibrated To Deliver; With Certificate of Traceability; White scale graduations; Capacity: 250mL; Subdivisions: 2mL	Each	4536
08557F	Fisherbrand* Cylinder; Serialized class A; With Certificate of Traceability; Borosilicate glass; Graduated; Calibrated to deliver; Capacity: 500mL; Subdivisions: 5mL; Tolerance: +/-2mL; Approx. height: 375mm	Each	5359
08557H	Fisherbrand* Cylinder; Serialized class A; With Certificate of Traceability; Borosilicate glass; Graduated; Calibrated to deliver; Capacity: 2000mL; Subdivisions: 20mL; Tolerance: +/-6mL; Approx. height: 520mm	Each	13330
0973030	Fisherbrand* Cylinder, Graduated; Serialized Class A Glass; Calibrated To Deliver; With Certificate of Traceability; White scale graduations; Capacity: 1000mL; Subdivisions: 10mL	Each	6502
08570A	Fisherbrand* Cylinders, Graduated; Class B Glass; Calibrated TD; Amber Scale; 5mL; Subdiv: 0.1mL	Each	1076
08570B	Fisherbrand* Cylinders, Graduated; Class B Glass; Calibrated TD; Amber Scale; 10mL; Subdiv: 0.2mL	Each	1154
08570C	Fisherbrand* Cylinders, Graduated; Class B Glass; Calibrated TD; Amber Scale; 25mL; Subdiv: 0.5mL	Each	1412
08570D	Fisherbrand* Cylinders, Graduated; Class B Glass; Calibrated TD; Amber Scale; 50mL; Subdiv: 1mL	Each	1466
08570E	Fisherbrand* Cylinders, Graduated; Class B Glass; Calibrated TD; Amber Scale; 100mL; Subdiv: 1mL	Each	1752
08570F	Fisherbrand* Cylinders, Graduated; Class B Glass; Calibrated TD; Amber Scale; 250mL; Subdiv: 2mL	Each	2421
08570G	Fisherbrand* Cylinders, Graduated; Class B Glass; Calibrated TD; Amber Scale; 500mL; Subdiv: 5mL	Each	3389
08570H	Fisherbrand* Cylinders, Graduated; Class B Glass; Calibrated TD; Amber Scale; 1000mL; Subdiv: 10mL	Each	5186

Part No.	Description	Pack	Price in Rs.
08570J	Fisherbrand* Cylinders, Graduated; Class B Glass; Calibrated TD; Amber Scale; 2000mL; Subdiv: 20mL	Each	8491
08550A	Fisherbrand* Cylinder; Standard class B; Borosilicate glass; Calibrated to deliver; White-ceramic graduations; Heavy round bases; Handy pouring spouts; Capacity: 5mL; Subdivisions: 0.1mL; Tolerance: +/-0.1mL; Approx. height: 120mm	Each	1100
08550B	Fisherbrand* Cylinder, Graduated; Class B; Borosilicate glass; Calibrated TD; White scale; Heavy round bases; Handy pouring spouts; Capacity: 10mL; Subdiv: 0.2mL; Tolerance: +/-0.20mL; Approx. height: 140mm	Each	1137
08550C	Fisherbrand* Cylinder, Graduated; Class B Glass; Calibrated To Deliver; White-ceramic graduations; Capacity: 25mL; Subdivisions: 0.5mL	Each	1398
08550D	Fisherbrand* Cylinder, Graduated; Class B Glass; Calibrated To Deliver; White-ceramic graduations; Capacity: 50mL; Subdivisions: 1mL	Each	1454
08550E	Fisherbrand* Cylinder, Graduated; Class B Glass; Calibrated To Deliver; White-ceramic graduations; Capacity: 100mL; Subdivisions: 1mL	Each	1734
08550F	Fisherbrand* Cylinder, Graduated; Class B; Borosilicate glass; Calibrated TD; White scale; Heavy round bases; Handy pouring spouts; Capacity: 250mL; Subdiv: 2mL; Tolerance: +/-2.0mL; Approx. height: 330mm	Each	2399
08550G	Fisherbrand* Cylinder, Graduated; Class B; Glass; Calibrated To Deliver; White-ceramic graduations; Handy pouring spouts; Capacity: 500mL; Subdivisions: 5mL; Tolerance: +/-4.0mL; Approx. height: 375mm	Each	3356
08550H	Fisherbrand* Cylinder, Graduated; Class B Glass; Calibrated To Deliver; White-ceramic graduations; Capacity: 1000mL; Subdivisions: 1mL	Each	5209
08550J	Fisherbrand* Cylinder, Graduated; Class B; Borosilicate glass; Calibrated TD; White scale; Heavy round bases; Handy pouring spouts; Capacity: 2000mL; Subdiv: 20mL; Tolerance: +/-12.0mL; Approx. height: 520mm	Each	8409
08555A	Fisherbrand* Cylinders, Graduated; Class B; Glass; With double ceramic metric scales; Calibrated To Deliver; Heavy base; Pouring lip; Capacity: 10mL; Subdivisions: 0.2mL; Tolerance: +/-0.1mL	Each	966
08555B	Fisherbrand* Cylinders, Graduated; Class B; Glass; With double ceramic metric scales; Calibrated To Deliver; Heavy base; Pouring lip; Capacity: 25mL; Subdivisions: 0.5mL; Tolerance: +/-0.20mL	Each	1023
08555C	Fisherbrand* Cylinders, Graduated; Class B; Glass; With double ceramic metric scales; Calibrated To Deliver; Heavy base; Pouring lip; Capacity: 50mL; Subdivisions: 1mL; Tolerance: +/-0.34mL	Each	1093
08555D	Fisherbrand* Cylinder, Graduated; Class B Glass; Calibrated To Deliver; Double metric graduations for reading up or down; Pouring lip; Capacity: 100mL; Subdivisions: 1mL; Tolerance: +/-1.0mL	Each	1224
08555E	Fisherbrand* Cylinder, Graduated; Class B; Borosilicate glass; Calibrated TD; Dual White Scale; Heavy base; Accuracy assured even at lowest graduations; Pouring lip; Capacity: 250mL; Subdiv: 2mL; Tolerance: +/-2.0mL; Case of: 18	Each	2303
08555F	Fisherbrand* Cylinder, Graduated; Class B Glass; Calibrated To Deliver; Double metric graduations for reading up or down; Pouring lip; Capacity: 500mL; Subdivisions: 5mL; Tolerance: +/-2.6mL	Each	2882
08555G	Fisherbrand* Cylinder, Graduated; Class B; Borosilicate glass; Calibrated TD; Dual White Scale; Heavy base; Accuracy assured even at lowest graduations; Pouring lip; Capacity: 1000mL; Subdiv: 10mL; Tolerance: +/-6.0mL; Case of: 8	Each	4745
S63455	Fisherbrand* Cylinders, Graduated; 10mL; 0.2mL Subdivisions; Borosilicate Glass; Single Metric Scale in Tough Ceramic Enamel; Heavy Glass Base; Rated To Deliver	Each	863
S63456	Fisherbrand* Cylinders, Graduated; 25mL; 0.5mL Subdivisions; Borosilicate Glass; Single Metric Scale in Tough Ceramic Enamel; Heavy Glass Base; Rated To Deliver	Each	779
S63457	Fisherbrand* Cylinders, Graduated; 50mL; 1mL Subdivisions; Borosilicate Glass; Single Metric Scale in Tough Ceramic Enamel; Heavy Glass Base; Rated To Deliver	Each	880
S63458	Fisherbrand* Cylinders, Graduated; 100mL; 1mL Subdivisions; Borosilicate Glass; Single Metric Scale in Tough Ceramic Enamel; Heavy Glass Base; Rated To Deliver	Each	1151
S63459	Fisherbrand* Cylinders, Graduated; 250mL; 2mL Subdivisions; Borosilicate Glass; Single Metric Scale in Tough Ceramic Enamel; Heavy Glass Base; Rated To Deliver	Each	1911
S63460	Fisherbrand* Cylinders, Graduated; 500mL; 5mL Subdivisions; Borosilicate Glass; Single Metric Scale in Tough Ceramic Enamel; Heavy Glass Base; Rated To Deliver	Each	2603
S63461	Fisherbrand* Cylinders, Graduated; 1000mL; 10mL Subdivisions; Borosilicate Glass; Single Metric Scale in Tough Ceramic Enamel; Heavy Glass Base; Rated To Deliver	Each	3769
08562A	Fisherbrand* Cylinder, Graduated; Glass; Low-form; Double metric scales for reading up or down; Capacity: 50mL; Subdivisions: 1mL	Each	1361
08562B	Fisherbrand* Cylinders, Graduated; Glass; Low-form; TD; Dual White Scale; Capacity: 100mL; Subdiv: 2mL	Each	1813

Part No.	Description	Pack	Price in Rs.
<b>Ungraduated Cylinders</b>			
08530A	Fisherbrand* Cylinder, Ungraduated; Borosilicate glass; With pouring lip; Approx. Capacity: 25mL; O.D.: 25mm	Each	915
08530B	Fisherbrand* Cylinder, Ungraduated; Borosilicate glass; With pouring lip; Approx. Capacity: 50mL; O.D.: 25mm	Each	973
08530C	Fisherbrand* Cylinder, Ungraduated; Borosilicate glass; With pouring lip; Approx. Capacity: 125mL; O.D.: 41mm	Each	1389
08530E	Fisherbrand* Cylinder, Ungraduated; Borosilicate glass; With pouring lip; Approx. Capacity: 175mL; O.D.: 41mm	Each	1483
08530F	Fisherbrand* Cylinder, Ungraduated; Borosilicate glass; With pouring lip; Approx. Capacity: 225mL; O.D.: 41mm	Each	1621
08530G	Fisherbrand* Cylinder, Ungraduated; Borosilicate glass; With pouring lip; Approx. Capacity: 275mL; O.D.: 41mm	Each	1675
08530J	Fisherbrand* Cylinder, Ungraduated; Borosilicate glass; With pouring lip; Approx. Capacity: 300mL; O.D.: 41mm	Each	1864
08530K	Fisherbrand* Cylinder, Ungraduated; Borosilicate glass; With pouring lip; Approx. Capacity: 350mL; O.D.: 41mm	Each	1897
08530L	Fisherbrand* Cylinder, Ungraduated; Borosilicate glass; With pouring lip; Approx. Capacity: 600mL; O.D.: 51mm	Each	2373
0973032	Fisherbrand* Cylinder, Ungraduated; Borosilicate glass; With pouring lip; Approx. Capacity: 1000mL; O.D.: 64mm	Each	3209
0973033	Fisherbrand* Cylinder, Ungraduated; Borosilicate glass; With pouring lip; Approx. Capacity: 2000mL; O.D.: 82.5mm	Each	6671
08535A	Fisherbrand* Cylinders, Ungraduated; Borosilicate glass; w/o pouring lip; 50mL; O.D.: 25mm	Each	1187
08535B	Fisherbrand* Cylinders, Ungraduated; Borosilicate glass; w/o pouring lip; 125mL; O.D.: 41mm	Each	1360
08535C	Fisherbrand* Cylinders, Ungraduated; Borosilicate glass; w/o pouring lip; 175mL; O.D.: 41mm	Each	1506
08535E	Fisherbrand* Cylinders, Ungraduated; Borosilicate glass; w/o pouring lip; 225mL; O.D.: 41mm	Each	2217
08535F	Fisherbrand* Cylinders, Ungraduated; Borosilicate glass; w/o pouring lip; 500mL; O.D.: 51mm	Each	2073
08535G	Fisherbrand* Cylinders, Ungraduated; Borosilicate glass; w/o pouring lip; 620mL; O.D.: 51mm	Each	2026
S40790	Fisherbrand* Cylinder, Ungraduated; 500mL; O.D. x H: 2 x 11.8 in.; Borosilicate glass; Heavy wall construction; Fused, round base increases stability	Each	1846
<b>Glass Graduated Cylinders for Special Purpose</b>			
FBA08558	Fisherbrand* Cylinders, Graduated; Glass; Petroleum Distillations; White scale; 100mL; Subdiv: 1mL	Each	1743
<b>Plastic Graduated Cylinders</b>			
115822	Fisherbrand* Cylinders, Hydrometer; Plastic; Features molded graduations, overflow cup, and 4.25 in. dia. (10cm) base; Capacity: 500mL; Subdivisions: 5mL	Each	1909
0300730	Fisherbrand* Cylinder, Graduated; Polymethylpentene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 10mL; Subdivisions: 0.2mL	Each	95
0300731	Fisherbrand* Cylinder, Graduated; Polymethylpentene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 25mL; Subdivisions: 0.5mL	Each	122
0300732	Fisherbrand* Cylinder, Graduated; Polymethylpentene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 50mL; Subdivisions: 1mL	Each	151
0300733	Fisherbrand* Cylinder, Graduated; Polymethylpentene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 100mL; Subdivisions: 1mL	Each	196
0300734	Fisherbrand* Cylinder, Graduated; Polymethylpentene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 250mL; Subdivisions: 2mL	Each	278
0300735	Fisherbrand* Cylinder, Graduated; Polymethylpentene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 500mL; Subdivisions: 5mL	Each	402
0300736	Fisherbrand* Cylinder, Graduated; Polymethylpentene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 1000mL; Subdivisions: 10mL	Each	680
0300737	Fisherbrand* Cylinder, Graduated; Polymethylpentene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 2000mL; Subdivisions: 20mL	Each	1957
0300738	Fisherbrand* Cylinder, Graduated; Polypropylene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 10mL; Subdivisions: 0.2mL	Each	105
0300739	Fisherbrand* Cylinder, Graduated; Polypropylene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 25mL; Subdivisions: 0.5mL	Each	113
0300740	Fisherbrand* Cylinder, Graduated; Polypropylene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 50mL; Subdivisions: 1mL	Each	132
0300741	Fisherbrand* Cylinder, Graduated; Polypropylene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 100mL; Subdivisions: 1mL	Each	157
0300742	Fisherbrand* Cylinder, Graduated; Polypropylene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 250mL; Subdivisions: 2mL	Each	229
0300743	Fisherbrand* Cylinder, Graduated; Polypropylene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 500mL; Subdivisions: 5mL	Each	295
0300744	Fisherbrand* Cylinder, Graduated; Polypropylene; 1000mL	Each	602
0300745	Fisherbrand* Cylinder, Graduated; Polypropylene; Calibrated To Contain/To Deliver; Autoclavable; Capacity: 2000mL; Subdivisions: 20mL	Each	1074