

## Digestion

Part No.	Description	Pack	Price in Rs.
FBA000188	Fisherbrand* Digestion System: Heat Shields for Tube Racks TC1000-4237, TC1000-0072	Each	9729

## Dishes

Part No.	Description	Pack	Price in Rs.
<b>Evaporating Dishes</b>			
02617139	Fisherbrand* Dishes, Evaporating; Tall Form; PTFE; 25mL; 38 Dia. x 34mmH; Withstands temperatures to 280 deg. C (536 deg. F)	Each	1403
02617140	Fisherbrand* Dishes, Evaporating; Tall Form; PTFE; 50mL; 40 Dia. x 50mmH; Withstands temperatures to 280 deg. C (536 deg. F)	Each	1738
02617141	Fisherbrand* Dishes, Evaporating; Tall Form; PTFE; 100mL; 64 Dia. x 50mmH; Withstands temperatures to 280 deg. C (536 deg. F)	Each	2541
02617142	Fisherbrand* Dishes, Evaporating; Tall Form; PTFE; 150mL; 78 Dia. x 43mmH; Withstands temperatures to 280 deg. C (536 deg. F)	Each	3303
02617144	Fisherbrand* Dishes, Evaporating; Tall Form; PTFE; 250mL; 97 Dia. x 50mmH; Withstands temperatures to 280 deg. C (536 deg. F)	Each	4041
02617145	Fisherbrand* Dishes, Evaporating; Tall Form; PTFE; 350mL; 100 Dia. x 60mmH; Withstands temperatures to 280 deg. C (536 deg. F)	Each	4490
02617146	Fisherbrand* Dishes, Evaporating; Low Form; PTFE; 25mL; 63 Dia. x 10mmH; Withstands temperatures to 280 deg. C (536 deg. F)	Each	1445
02617147	Fisherbrand* Dishes, Evaporating; Low Form; PTFE; 50mL; 79 Dia. x 13mmH; Withstands temperatures to 280 deg. C (536 deg. F)	Each	1831
02617148	Fisherbrand* Dishes, Evaporating; Low Form; PTFE; 100mL; 100 Dia. x 16mmH; Withstands temperatures to 280 deg. C (536 deg. F)	Each	1375
02617149	Fisherbrand* Dishes, Evaporating; Low Form; PTFE; 180mL; 121 Dia. x 19mmH; Withstands temperatures to 280 deg. C (536 deg. F)	Each	3620
02617150	Fisherbrand* Dishes, Evaporating; Low Form; PTFE; 350mL; 130 Dia. x 34mmH; Withstands temperatures to 280 deg. C (536 deg. F)	Each	4820
02617151	Fisherbrand* Dishes, Evaporating; Low Form; PTFE; 400mL; 130 Dia. x 36mmH; Withstands temperatures to 280 deg. C (536 deg. F)	Each	3517
<b>Weighing Dishes</b>			
08732105	Fisherbrand* Dishes, Weighing; Aluminum; Disposable; 20mL; Slightly Fluted Sides with Flanged Rims	Pack of 100	382
08732106	Fisherbrand* Dishes, Weighing; Aluminum; Disposable; 60mL; Slightly Fluted Sides with Flanged Rims	Pack of 100	486
08732107	Fisherbrand* Dishes, Weighing; Aluminum; Disposable; 75mL; Slightly Fluted Sides with Flanged Rims	Pack of 100	434
08732108	Fisherbrand* Dishes, Weighing; Aluminum; Disposable; 200mL; Slightly Fluted Sides with Flanged Rims	Pack of 50	472
08732109	Fisherbrand* Dishes, Weighing; Aluminum; Disposable; 500mL; Slightly Fluted Sides with Flanged Rims	Pack of 50	658
08732110	Fisherbrand* Dishes, Evaporating; Aluminum Weighing/Drying Pan, Disposable; Top ID: 12.2cm, with flat top	Pack of 100	874
08732111	Fisherbrand* Dishes, Evaporating; Aluminum Weighing/Drying Pan, Disposable; Top ID: 10.9cm, w/fluted sides	Pack of 50	654
<b>Glass Special Purpose Dishes</b>			
FBA08807	Fisherbrand* Dishes, Special Purpose; Preparation; Stackable; Glass; 300mL; Dia. x H: 4-1/2 x 2 in.	Each	608
<b>Plastic Petri Dishes</b>			
08757150	Fisherbrand* Petri Dish, Disposable; Sectioned plates; W/numbered sections; Radiation sterilized; Polystyrene; 25/polyethylene sleeve; Size: 100 dia. x 15mm H; Type: I (2 sections)	Case of 500	6866
087582	Fisherbrand* Petri Dish, Disposable; Sectioned plates; W/numbered sections; Radiation sterilized; Polystyrene; 25/polyethylene sleeve; Size: 100 dia. x 15mm H; Type: X (4 sections)	Case of 500	7281
087583	Fisherbrand* Petri Dish, Disposable; Sectioned plates; W/numbered sections; Radiation sterilized; Polystyrene; 25/polyethylene sleeve; Size: 100 dia. x 15mm H; Type: Y (3 sections)	Case of 500	7201
0875711C	Fisherbrand* Plates, Contact; Disposable; 65 x 15mm Sterile; Sixteen 10 x 10mm square grid; Marked numerically and alphabetically	Case of 500	6179
09720500	Fisherbrand* Dish, Disposable petri; Molded polystyrene; Dia.: 47mm; Sterilized; Dishes only; Tight fitting lid prevents dehydration; 47mm dishes can be used with solid agar medium;	Pack of 150	3681
09720501	Fisherbrand* Dish, Disposable petri; Molded polystyrene; Dia.: 47mm; Sterilized; Dishes w/pads; Minimizes contamination; Pure cellulose pad; Extended 3-year shelf-life;	Pack of 150	4801
09720502	Fisherbrand* Dish, Disposable petri; Molded polystyrene; Dia.: 47mm; Sterilized; Dishes only; Tight fitting lid prevents dehydration; 47mm dishes can be used with solid agar medium;	Pack of 600	10721

Part No.	Description	Pack	Price in Rs.
09720503	Fisherbrand* Dish, Disposable petri; Molded polystyrene; Dia. 47mm; Sterilized; Dishes w/pads; Minimizes contamination; Pure cellulose pad; Extended 3-year shelf-life;	Pack of 600	14722
0875711	Fisherbrand* Petri Dish, Disposable; Polystyrene; Vented lid; Extra-deep for deeper media fills; Ideal for fungal cultures, plant propagation; radiation sterilized; 20 per polyethylene sleeve; Dia. x H: 100 x 25mm	Case of 500	14952
0875711Z	Fisherbrand* Petri Dish, Disposable; Polystyrene; Stackable; Extra-deep for deeper media fills; Ideal for fungal cultures, plant propagation; radiation sterilized; 20 per polyethylene sleeve; Dia. x H: 100 x 20mm	Case of 500	12524
0875711YZ	Fisherbrand* Petri Dish, Disposable; Media-Miser; Polystyrene; Smaller size saves media; Clear tops; Snug fit lids; Venting ribs; Raised, straight ridge for stacking; Radiation-sterilized; Dia. x H (approx): 35 x 10mm	Case of 500	10545
0875713A	Fisherbrand* Petri Dish, Disposable; Media-Miser; Polystyrene; Smaller size saves media; Clear tops; Snug fit lids; Venting ribs; Raised, straight ridge for stacking; Radiation-sterilized; Dia. x H (approx): 60 x 15mm	Case of 500	4993
0875714G	Fisherbrand* Petri Dish, Disposable; Media-Miser; Polystyrene; Smaller size saves media; Clear tops; Snug fit lids; Venting ribs; Raised, straight ridge for stacking; Radiation-sterilized; Dia. x H (approx): 95 x 15mm	Case of 500	6220
S67961	Fisherbrand* Petri Dish, Disposable; Media-Miser; Polystyrene; Smaller size saves media; Clear tops; Snug fit lids; Venting ribs; Raised, straight ridge for stacking; Radiation-sterilized; Dia. x H (approx): 60 x 15mm	Pack of 20	200
S67962	Fisherbrand* Media-Miser Petri Dishes FSE Diameter x Height (approximate): 85 x 15mm; Percent Media Saved: 11%	Pack of 25	312
S68051	Fisherbrand* Petri Dish, Partitioned; Disposable; Stackable; One hand operation; Ventilation ribs; Radiation sterilized; Quad-plate; X-type; 100 X 15mm;	Pack of 25	363
S68052	Fisherbrand* Petri Dish, Partitioned; Disposable; Stackable; One hand operation; Ventilation ribs; Radiation sterilized; Tri-plate; Y-Type; 100 X 15mm;	Pack of 25	361
S68053	Fisherbrand* Petri Dish, Partitioned; Disposable; Stackable; One hand operation; Ventilation ribs; Radiation sterilized; Bi-plate; I-Type; 100 X 15mm;	Pack of 25	344
0875712	Fisherbrand* Dish, Petri; Clear Lid; Radiation-sterilized; Polystyrene; Raised Ridge; 100 Dia. x 15mmH; Edges are square; Three venting ribs in lid; Stackable; Disposable	Case of 500	5759
0875713	Fisherbrand* Dish, Petri; Clear Lid; Radiation-sterilized; Polystyrene; Beveled Ridge; 100 O.D. x 15mm H; Edges are square; Three venting ribs in lid; Stackable; Disposable	Case of 500	5759
0875714	Fisherbrand* Dish, Petri; Clear Lid; Radiation-sterilized; Polystyrene; Beveled Ridge; 150 Dia. x 15mmH; Edges are square; Three venting ribs in lid; Stackable; Disposable	Case of 100	4645
087579B	Fisherbrand* Dish, Petri; Disposable; Stackable; Vented; Radiation sterilized; Raised Ridge; 100 x 15mm	Case of 500	6754
0875711A	Fisherbrand* Petri Dish, Square; Disposable; W/Grid of thirty-six 13 x 13mm squares; Polystyrene; Radiation Sterilized; L x W x H: 100 x 100 x 15mm; 10 per polyethylene sleeve	Case of 500	16801
<b>Petri Dish Racks</b>			
1479300	Fisherbrand* Dishes, Petri, Accessories; Rack; Easyplate; 9 in. (22.8cm) rack holds 48 100mm petri dishes; Color: Natural; Polycarbonate; Autoclavable; Chemical-Resistant	Each	12888
1479301	Fisherbrand* Dishes, Petri, Accessories; Rack; Easyplate; 9 in. (22.8cm) rack holds 48 100mm petri dishes; Color: Red; Polycarbonate; Autoclavable; Chemical-Resistant	Each	12888
1479302	Fisherbrand* Dishes, Petri, Accessories; Rack; Easyplate; 9 in. (22.8cm) rack holds 48 100mm petri dishes; Color: Yellow; Polycarbonate; Autoclavable; Chemical-Resistant	Each	12888
1479303	Fisherbrand* Dishes, Petri, Accessories; Rack; Easyplate; 9 in. (22.8cm) rack holds 48 100mm petri dishes; Color: Blue; Polycarbonate; Autoclavable; Chemical-Resistant	Each	12888
1479304	Fisherbrand* Dishes, Petri, Accessories; Rack; Easyplate; 9 in. (22.8cm) rack holds 48 100mm petri dishes; Colors: Natural, Red, Yellow, Blue; Polycarbonate; Autoclavable; Chemical-Resistant	Pack of 4	49468
1479305	Fisherbrand* Dishes, Petri, Accessories; Rack; Easyplate; 14 in. (35.5cm) rack holds 84 100mm petri dishes; Colors: Natural; Polycarbonate; Autoclavable; Chemical-Resistant	Each	14375
1479306	Fisherbrand* Dishes, Petri, Accessories; Rack; Easyplate; 14 in. (35.5cm) rack holds 84 100mm petri dishes; Colors: Blue; Polycarbonate; Autoclavable; Chemical-Resistant	Each	14375
1479317	Fisherbrand* Dishes, Petri, Accessories; Rack for Rodak plates; Clear acrylic; Holds twenty-eight 65mm dia. Plates; Contains solid acrylic, permanently fixed posts to prevent shifting of plates and slide-down cover	Each	3596
<b>Accessories for Petri Dishes</b>			
1479307	Fisherbrand* Dishes, Petri, Accessories; Incubation Tray; Stackable; Hold five 100mm or fifteen 60mm Petri dishes; Chemical-resistant polypropylene; Dimensions: 9.88L x 9.31W x 1.38 in. H (25 x 24 x 3cm)	Pack of 3	2994
<b>Boxes</b>			
FBA03460	Fisherbrand* Dishes, Petri, Accessories; Stainless Steel Sterilizer Box; Holds 18 Dishes	Each	2240