

Electrical Equipment

Part No.	Description	Pack	Price in Rs.
Electrical Equipment			
012411	Fisherbrand* Multimeter, Digital Auto-Range, Ranges from 1 millivolt to 400 volts AC/DC, Audible beep for checking continuity, with case, batteries, probes and instructions	Each	3293
1464816	Fisherbrand* Traceable Battery Tester, checks voltage 1.35 to 12V, Tests in one second, Self-powered, No internal battery	Each	2532
Electrical Supplies			
095205	Fisherbrand* Surge Suppressor Strip	Each	4535

Electroporators

Part No.	Description	Pack	Price in Rs.
Cuvets for Electroporation			
FBAFB101	Fisherbrand* Cuvettes, Electroporation; Cuvettes Plus; 1mm gap with imbedded aluminum electrodes, Capacity: 90uL; Field strengths up to 25 kV/cm, gamma irradiated for guaranteed sterility, Disposable	Pack of 50	8916
FBAFB102	Fisherbrand* Cuvettes, Electroporation Cuvettes Plus; 2mm gap with imbedded aluminum electrodes, Capacity: 400microL; Field strengths up to 25 kV/cm, Gamma irradiated for guaranteed sterility, Disposable	Pack of 50	8916
FBAFB104	Fisherbrand* Cuvettes, Electroporation Cuvettes Plus; 4mm gap with imbedded aluminum electrodes, Capacity: 800uL; Field strengths up to 25 kV/cm, gamma irradiated for guaranteed sterility, Disposable	Pack of 50	8916

Filtration

Part No.	Description	Pack	Price in Rs.
Qualitative Filter Paper			
FB59043	Fisherbrand* Filter paper QL 110 qualitative medium/fine retention and medium flow rate cellulose equivalent to Whatman No 2 125mm diameter	Pack of 100	699
FB59045	Fisherbrand* Filter paper QL 110 qualitative medium/fine retention and medium flow rate cellulose equivalent to Whatman No 2 150mm diameter	Pack of 100	820
FB59047	Fisherbrand* Filter paper QL 110 qualitative medium/fine retention and medium flow rate cellulose equivalent to Whatman No 2 185mm diameter	Pack of 100	1069
FB59049	Fisherbrand* Filter paper QL 110 qualitative medium/fine retention and medium flow rate cellulose equivalent to Whatman No 2 240mm diameter	Pack of 100	1704
FB59051	Fisherbrand* Filter paper QL 110 qualitative medium/fine retention and medium flow rate cellulose equivalent to Whatman No 2 270mm diameter	Pack of 100	2603
FB59055	Fisherbrand* Filter paper QL 110 qualitative medium/fine retention and medium flow rate cellulose equivalent to Whatman No 2 460mm x 570mm 100 sheets	Pack of 100	7002
FB59041	Fisherbrand* Filter paper QL 110 qualitative medium/fine retention and medium flow rate cellulose equivalent to Whatman No 2 70mm diameter	Pack of 100	380
FB59061	Fisherbrand* Filter paper QL 115 qualitative same as QL110 already folded cellulose equivalent to Whatman No 2V 125mm diameter	Pack of 100	2169
FB59063	Fisherbrand* Filter paper QL 115 qualitative same as QL110 already folded cellulose equivalent to Whatman No 2V 150mm diameter	Pack of 100	2324
FB59065	Fisherbrand* Filter paper QL 115 qualitative same as QL110 already folded cellulose equivalent to Whatman No 2V 185mm diameter	Pack of 100	2573
FB59067	Fisherbrand* Filter paper QL 115 qualitative same as QL110 already folded cellulose equivalent to Whatman No 2V 240mm diameter	Pack of 100	3594
FB59081	Fisherbrand* Filter paper QL 120 qualitative medium/fine retention and high loading capacity cellulose equivalent to Whatman No 3 125mm diameter	Pack of 100	1115
FB59083	Fisherbrand* Filter paper QL 120 qualitative medium/fine retention and high loading capacity cellulose equivalent to Whatman No 3 150mm diameter	Pack of 100	2093
FB59075	Fisherbrand* Filter paper QL 120 qualitative medium/fine retention and high loading capacity cellulose equivalent to Whatman No 3 70mm diameter	Pack of 100	680
FB59077	Fisherbrand* Filter paper QL 120 qualitative medium/fine retention and high loading capacity cellulose equivalent to Whatman No 3 90mm diameter	Pack of 100	854
FB59099	Fisherbrand* Filter paper QL 125 qualitative coarse retention and fast flow rate cellulose equivalent to Whatman No 4 110mm diameter	Pack of 100	559
FB59101	Fisherbrand* Filter paper QL 125 qualitative coarse retention and fast flow rate cellulose equivalent to Whatman No 4 125mm diameter	Pack of 100	699
FB59103	Fisherbrand* Filter paper QL 125 qualitative coarse retention and fast flow rate cellulose equivalent to Whatman No 4 150mm diameter	Pack of 100	820