

Surface Tension

Part No.	Description	Price in Rs.
FBA14818	Fisher Scientific* Surface Tension Apparatus	8357
148125	Fisher Scientific* Test Ring; Tensiometer Replacement Platinum-Iridium; Test ring	34301
14818A	Fisher Scientific* Surface Tension Apparatus Replacement Capillary Tubes	6000

Ultrasonic Cleaners

Part No.	Description	Price in Rs.
FB11002	Fisherbrand* Ultrasonic bath digital, 2.75L capacity	POR
FB11139	Fisherbrand* Cover for ultrasonic bath for FB11002 stainless steel	4557
FB11144	Fisherbrand* Basket for ultrasonic bath for FB11002 stainless steel	7583
FB11003	Fisherbrand* Ultrasonic bath digital, 6.9L capacity	POR
FB11140	Fisherbrand* Cover for ultrasonic bath for FB11003 stainless steel	9861
FB11145	Fisherbrand* Basket for ultrasonic bath for FB11003 stainless steel	25026
FB11004	Fisherbrand* Ultrasonic bath digital, 9.50L capacity	POR
FB11141	Fisherbrand* Cover for ultrasonic bath for FB11004 / FB11005 stainless steel	6817
FB11146	Fisherbrand* Basket for ultrasonic bath for FB11004 stainless steel	22930
FB11005	Fisherbrand* Ultrasonic bath digital, 12.75L capacity	POR
FB11147	Fisherbrand* Basket for ultrasonic bath for FB11005 stainless steel	22930
FB11006	Fisherbrand* Ultrasonic bath digital, 18.00L capacity	POR
FB11142	Fisherbrand* Cover for ultrasonic bath for FB11006 stainless steel	7965
FB11148	Fisherbrand* Basket for ultrasonic bath for FB11006 stainless steel	25974
FB11007	Fisherbrand* Ultrasonic bath digital, 28.00L capacity	POR
FB11143	Fisherbrand* Cover for ultrasonic bath for FB11007 stainless steel	10809
FB11149	Fisherbrand* Basket for ultrasonic bath for FB11007 stainless steel	29765
FB11052	Fisherbrand* Cover for ultrasonic bath for FB11028 stainless steel	11848
FB11082	Fisherbrand* Basket ultrasonic bath for FB11074 stainless steel 100mm x 456mm x 46mm	25026
FB11087	Fisherbrand* Cover ultrasonic bath stainless steel for FB11074	9861
FB11172	Fisherbrand* Bath TI-H-25 MF3 dual frequency 230V 50Hz 2300W 19.8L 35/130 kHz	POR
FB11173	Fisherbrand* Bath TI-H-25 MF2 dual frequency 230V 50Hz 2300W 19.8L 25/45 kHz	POR
FB11175	Fisherbrand* Bath TI-H-55 MF3 dual frequency 400V 50Hz 3600W 3 phase 45.0L 35/130 kHz	POR
FB11178	Fisherbrand* Bath TI-H-80 MF3 dual frequency 400V 50Hz 4900W 3 phase 67.5L 35/130 kHz	POR
FB11179	Fisherbrand* Bath TI-H-80 MF2 dual frequency 400V 50Hz 4900W 3 phase 67.5L 25/45 kHz	POR
FB11181	Fisherbrand* Bath TI-H-115 MF3 dual frequency 400V 50Hz 7400W 3 phase 90.0L 35/130 kHz	POR
FB11182	Fisherbrand* Bath TI-H-115 MF dual frequency 400V 50Hz 7400W 3 phase 90.0L 25/45 kHz	POR
FB11184	Fisherbrand* Bath TI-H-160 MF3 dual frequency 400V 50Hz 7800W 3 phase 135L 35/130 kHz	POR
FB11185	Fisherbrand* Bath TI-H-160 MF2 dual frequency 400V 50Hz 7800W 3 phase 135L 25/45 kHz	POR
FB11035	Fisherbrand* Cover for ultrasonic bath for FB15046/FB15047 stainless steel	4375
FB15070	Fisherbrand* Basket ultrasonic bath for FB15046/15047	5833
FB15081	Fisherbrand* Cover ultrasonic bath for FB15046/15047	2187
FB15071	Fisherbrand* Basket ultrasonic bath for FB15048/FB15049	6380

Prices mentioned are for one unit unless pack size is mentioned in the description

Contact: India Toll Free.: 1 800 209 7001 Fax: 022-6680 3001 / 02 e-mail: qfc.customer@thermofisher.com

Part No.	Description	Price in Rs.
FB15082	Fisherbrand* Cover ultrasonic bath for FB15048/FB15049	3099
FB11037	Fisherbrand* Cover for ultrasonic bath for FB15050/FB15051/FB15052/FB15053 stainless steel	7017
FB11041	Fisherbrand* Cover for ultrasonic bath for FB15050/FB15051/FB15052/FB15053 transparent	5286
FB15083	Fisherbrand* Cover ultrasonic bath for FB15050/FB15051/FB15052/FB15053	3463
FB15072	Fisherbrand* Basket ultrasonic bath for FB15050/FB15051	7473
FB11037	Fisherbrand* Cover for ultrasonic bath for FB15050/FB15051/FB15052/FB15053 stainless steel	7017
FB11041	Fisherbrand* Cover for ultrasonic bath for FB15050/FB15051/FB15052/FB15053 transparent	5286
FB15083	Fisherbrand* Cover ultrasonic bath for FB15050/FB15051/FB15052/FB15053	3463
FB11037	Fisherbrand* Cover for ultrasonic bath for FB15050/FB15051/FB15052/FB15053 stainless steel	7017
FB11041	Fisherbrand* Cover for ultrasonic bath for FB15050/FB15051/FB15052/FB15053 transparent	5286
FB15083	Fisherbrand* Cover ultrasonic bath for FB15050/FB15051/FB15052/FB15053	3463
FB15073	Fisherbrand* Basket ultrasonic bath for FB15052/FB15053	8567
FB15074	Fisherbrand* Basket ultrasonic bath for FB15054/FB15055	10754
FB15084	Fisherbrand* Cover ultrasonic bath for FB15054/FB15055	4375
FB15075	Fisherbrand* Basket stainless steel 465mm x 106mm x 50 mm for use with FB15056/FB15057	24971
FB15076	Fisherbrand* Basket stainless steel 455mm x 106mm x 75 mm for use with FB15058/FB15059	29346
FB11060	Fisherbrand* Basket for ultrasonic bath for FB15060/FB15061 stainless steel	22930
FB11055	Fisherbrand* Cover for ultrasonic bath for FB15060/FB15061/FB15062/FB15063 transparent	8020
FB15077	Fisherbrand* Basket stainless steel 300mm x 240mm x 150mm for use with FB15060/FB15061	22784
FB11060	Fisherbrand* Basket for ultrasonic bath for FB15060/FB15061 stainless steel	22930
FB15086	Fisherbrand* Cover ultrasonic bath for FB15060/FB15061/FB15062/FB15063	5468
FB15086	Fisherbrand* Cover ultrasonic bath for FB15060/FB15061/FB15062/FB15063	5468
FB11055	Fisherbrand* Cover for ultrasonic bath for FB15060/FB15061/FB15062/FB15063 transparent	8020
FB11061	Fisherbrand* Basket ultrasonic bath for FB15062/FB15063	22930
FB11055	Fisherbrand* Cover for ultrasonic bath for FB15060/FB15061/FB15062/FB15063 transparent	8020
FB15086	Fisherbrand* Cover ultrasonic bath for FB15060/FB15061/FB15062/FB15063	5468
FB15065	Fisherbrand* Ultrasonic bath heated 230V 50Hz 18.00L	POR
FB15079	Fisherbrand* Basket ultrasonic bath for FB15064/FB15065	25883
FB15087	Fisherbrand* Cover ultrasonic bath for FB15064/FB15065	5833
FB15080	Fisherbrand* Basket ultrasonic bath for FB15066/FB15067	29710
FB15088	Fisherbrand* Cover ultrasonic bath for FB15066/FB15067	6926
FB15069	Fisherbrand* Bath S900H ultrasonic heated 715mm (w) x 570mm (d) x 467mm (h) 230V 50/60Hz 2800W 90L	POR
15335103	Fisher Scientific* Cleaner, Ultrasonic Digital Control; FS110D; With heater; Maintains 60 deg. C temperature; 10 qt. (9.5L) capacity; 120V 50/60Hz, 4A, 469w; Exterior L x W x H: 12.25 x 12 x 11.25 in. (31.1 x 30.5 x 28.6cm)	POR
15335104	Fisher Scientific* Cleaner, Ultrasonic Digital Control; FS140D; With heater; Maintains 60 deg. C temperature; 13 qt. (12.3L) capacity; 120V 50/60Hz, 4A, 510w; Exterior L x W x H: 12.25 x 12.25 x 12.75 in. (31.1 x 31.1 x 32.4cm)	POR
15335105	Fisher Scientific* Cleaner, Ultrasonic Digital Control; FS220D; With heater; Maintains 60 deg. C temperature; 20 qt. (18.9L) capacity; 120V 50/60Hz, 8A, 881w; Exterior L x W x H: 14.75 x 20.25 x 11.25 in. (37.5 x 52.1 x 28.6cm)	POR
15335100	Fisher Scientific* Cleaner, Ultrasonic Digital Control; FS20D; With heater; Maintains 60 deg C temperature; 3 qt. (2.8L) capacity; 120V 50/60Hz, 1A, 143w; Exterior L x W x H: 8 x 9 x 11.25 in. (20.3 x 22.9 x 28.6cm)	79203
15335101	Fisher Scientific* Cleaner, Ultrasonic Digital Control; FS30D; With heater; Maintains 60 deg. C temperature; 4 qt. (3.8L) capacity; 120V 50/60Hz, 2A, 239w; Exterior L x W x H: 8.5 x 10 x 11.25 in. (21.6 x 25.4 x 28.6cm)	POR
15335102	Fisher Scientific* Cleaner, Ultrasonic Digital Control; FS60D; With heater; Maintains 60 deg. C temperature; 6 qt. (5.7L) capacity; 120V 50/60Hz, 3A; Exterior L x W x H: 8.5 X2.25X11.25 in. (21.6 x 31.1 x 28.6cm)	POR
15335110	Fisher Scientific* Cleaner, Ultrasonic 10qt.Cap (9.5 L) Capacity; 120V 50/60Hz, 2A, 185w; Without heater; Stainless-steel cover; Overall (Exterior) L x W x H: 12.25 x 12 x 11.25 in. (31.1 x 30.5 x 28.6cm)	POR

Prices mentioned are for one unit unless pack size is mentioned in the description

Part No.	Description	Price in Rs.
15335112	Fisher Scientific* Cleaner, Ultrasonic 10qt. (9.5 L)Capacity; 120V 50/60Hz, 4A, 469w; With heater; Stainless-steel cover; Overall (Exterior) L x W x H: 9.5 x 11.75 x 6 in. (24.1 x 29.8 x 15.2cm)	POR
15335140	Fisher Scientific* Cleaners, Ultrasonic 13 qt. (12.3L) Capacity; 120V 50/60Hz, 2A, 210w; Overall (Exterior) L x W x H: 12.25 x 12.25 x 12.75 in. (31.1 x 31.1 x 32.4cm); Without heater	POR
153356	Fisher Scientific* Cleaners, Ultrasonic Without heater; 18 oz.(0.5L) Capacity; 120V 50/60Hz, 0.4A, 40w	23986
1533520	Fisher Scientific* Cleaners, Ultrasonic Without heater; 3 qt. (2.8L) capacity; 60 min. timer; 120V 50/60Hz, 1A; 80w; 8 L x 9 W x 11-1/4 in. H; Full power instantly; No warmup; Stainless-steel tank and cover	52564
1533560	Fisher Scientific* Cleaners, Ultrasonic Without heater; 6 qt. (5.7L) capacity; 60 min. timer; 120V 50/60Hz, 1A; 150w; 8.5L x 12.25 W x 11.25 H in. (21.6 x 31.1 x 28.6cm); Stainless-steel tank and cover	96258
1533562	Fisher Scientific* Cleaners, Ultrasonic With heater; 6 qt. (5.7L) Capacity; 120V 50/60Hz, 3A; 335w; 8.5 L x 12.25 W x 11.25 in. H; Full power instantly, No warmup; Stainless-steel tank and cover	POR
15335220	Fisher Scientific* Cleaners, Ultrasonic Without heater; 22 qt. (18.9L) Capacity; 120V 50/60Hz, 3A, 320w; 14.75L x 20.5 Wx 11.25 H in. (37.5 x 52.1 x 28.6cm)	POR
1533530	Fisher Scientific* Cleaners, Ultrasonic Without heater; 4 qt. (3.8L) capacity; 60 min. timer; 120V 50/60Hz, 1A; 130w; 8.5 L x 10 Wx 11.25 H in. (21.6 x 25.4 x 28.6cm); Stainless-steel tank and cover	68160
15335222	Fisher Scientific* Cleaners, Ultrasonic With heater; 22 qt. (18.9L) Capacity; 120V 50/60Hz, 8A, 881w; 14.75L x 20.5W x 11.25 H in. (37.5 x 52.1 x 28.6cm)	POR
1533522	Fisher Scientific* Cleaners, Ultrasonic With heater; 3 qt (2.8L) Capacity; 120V 50/60Hz, 1A, 143w; 8Lx 9 W x 11.25 H in. (20.3 x 22.9 x 28.6cm)	57987
1533532	Fisher Scientific* Cleaners, Ultrasonic; ; With heater; Maintains 60 deg C temperature; 4 qt.(3.8L) capacity, 60 min. timer; 120V 50/60Hz, 2A, 239w; 8-1/2 L x 10 W x 11-1/4 in. H; Full power instantly, No warmup; Stainless-steel tank and cover	79041
15335142	Fisher Scientific* Cleaners, Ultrasonic 13 qt. (12.3L) Capacity; 120V 50/60Hz, 4A, 510w; Overall (Exterior) L x W x H: 12.25 x 12.25 x 12.75 in. (31.1 x 31.1 x 32.4cm); With heater	POR
LC0015	Fisher Scientific* Qualigens* ULTRASONIC CLEANER 1.5 lit	17000
LC0035	Fisher Scientific* Qualigens* ULTRASONIC CLEANER 3.5 lit	27000
LC0065	Fisher Scientific* Qualigens* ULTRASONIC CLEANER 6.5 lit	39000
LC0090	Fisher Scientific* Qualigens* ULTRASONIC CLEANER 9.0 lit	47000
FB15046	Fisherbrand* Ultrasonic bath,Elmasonic S 10/(Non Heated) Ultrasonic frequency: 37 kHz Tank capacity: 0,8 l Tank internal dimensions (W/D/H): 190x85x60 mm for use with S 10	27435
FB15047	Fisherbrand* Ultrasonic bath,Elmasonic S 10/(H) Ultrasonic frequency: 37 kHz Tank capacity: 0,8 l Tank internal dimensions (W/D/H): 190x85x60 mm for use with S 10 H	33481
FB15048	Fisherbrand* Ultrasonic bath,Elmasonic S 15/ Ultrasonic frequency: 37 kHz Tank capacity: 1,75 l Tank internal dimensions (W/D/H): 151x137x100 mm for use with S 15	49270
FB15049	Fisherbrand* Ultrasonic bath,Elmasonic S 15/(H) Ultrasonic frequency: 37 kHz Tank capacity: 1,75 l Tank internal dimensions (W/D/H): 151x137x100 mm for use with S15 H	58788
FB15050	Fisherbrand* Ultrasonic bath,Elmasonic S 30/ Ultrasonic frequency: 37 kHz Tank capacity: 2,75 l Tank internal dimensions (W/D/H): 240x137x100 mm for use with 30	55877
FB15051	Fisherbrand* Ultrasonic bath,Elmasonic S 30/(H) Ultrasonic frequency: 37 kHz Tank capacity: 2,75 l Tank internal dimensions (W/D/H): 240x137x100 mm for use with S 30H	58788
FB15052	Fisherbrand* Ultrasonic bath,Elmasonic S 40/ Ultrasonic frequency: 37 kHz Tank capacity: 4,25 l Tank internal dimensions (W/D/H): 240x137x150 mm for use with S 40	80176
FB15053	Fisherbrand* Ultrasonic bath,Elmasonic S 40/(H) Ultrasonic frequency: 37 kHz Tank capacity: 4,25 l Tank internal dimensions (W/D/H): 240x137x150 mm for use with S 40H	93166
FB15054	Fisherbrand* Ultrasonic bath,Elmasonic S 60/ Ultrasonic frequency: 37 kHz Tank capacity: 5,75 l Tank internal dimensions (W/D/H): 300x151x150 mm for use with S 60	92942
FB15055	Fisherbrand* Ultrasonic bath,Elmasonic S 60/(H) Ultrasonic frequency: 37 kHz Tank capacity: 5,75 l Tank internal dimensions (W/D/H): 300x151x150 mm for use with S 60 H	POR
FB15056	Fisherbrand* Ultrasonic bath,Elmasonic S 70/ Ultrasonic frequency: 37 kHz Tank capacity: 6,9 l Tank internal dimensions (W/D/H): 505x137x100 mm for use with S 70	POR
FB15057	Fisherbrand* Ultrasonic bath,Elmasonic S 70/(H) Ultrasonic frequency: 37 kHz Tank capacity: 6,9 l Tank internal dimensions (W/D/H): 505x137x100 mm for use with S 70 H	POR
FB15058	Fisherbrand* Ultrasonic bath,Elmasonic S 80/ Ultrasonic frequency: 37 kHz Tank capacity: 9,4 l Tank internal dimensions (W/D/H): 505x137x150 mm.Product profile S 80 / H for use with S 80	POR

Prices mentioned are for one unit unless pack size is mentioned in the description

Contact: India Toll Free.: 1 800 209 7001 Fax: 022-6680 3001 / 02 e-mail: qfc.customer@thermofisher.com

Part No.	Description	Price in Rs.
FB15059	Fisherbrand* Ultrasonic bath, Elmasonic S 80/(H) Ultrasonic frequency: 37 kHz Tank capacity: 9,4 l Tank internal dimensions (W/D/H): 505x137x150 mm. Product profile S 80 / H for use with S 80 H	POR
FB15060	Fisherbrand* Ultrasonic bath, Elmasonic S 100/ Ultrasonic frequency: 37 kHz Tank capacity: 9,5 l Tank internal dimensions (W/D/H): 300x240x150 mm for use with S 100	POR
FB15061	Fisherbrand* Ultrasonic bath, Elmasonic S 100/(H) Ultrasonic frequency: 37 kHz Tank capacity: 9,5 l Tank internal dimensions (W/D/H): 300x240x150 mm for use with S100 H	POR
FB15062	Fisherbrand* Ultrasonic bath, Elmasonic S 120/ Ultrasonic frequency: 37 kHz Tank capacity: 12,75 l Tank internal dimensions (W/D/H): 300x240x200 mm for use with S 120	POR
FB15063	Fisherbrand* Ultrasonic bath, Elmasonic S 120/(H) Ultrasonic frequency: 37 kHz Tank capacity: 12,75 l Tank internal dimensions (W/D/H): 300x240x200 mm for use with S120 H	POR
FB15064	Fisherbrand* Ultrasonic bath, Elmasonic S 180/ Ultrasonic frequency: 37 kHz Tank capacity: 18 l Tank internal dimensions (W/D/H): 327x300x200 mm for use with S 180	POR
FB15065	Fisherbrand* Ultrasonic bath, Elmasonic S 180/(H) Ultrasonic frequency: 37 kHz Tank capacity: 18 l Tank internal dimensions (W/D/H): 327x300x200 mm for use with S 180 H	POR
FB15066	Fisherbrand* Ultrasonic bath, Elmasonic S 300/ Ultrasonic frequency: 37 kHz Tank capacity: 28 l Tank internal dimensions (W/D/H): 505x300x200 mm for use with S 300	POR
FB15067	Fisherbrand* Ultrasonic bath, Elmasonic S 300/(H) Ultrasonic frequency: 37 kHz Tank capacity: 28 l Tank internal dimensions (W/D/H): 505x300x200 mm for use with S 300 H	POR
FB15068	Fisherbrand* Ultrasonic bath, Elmasonic S 450/ Ultrasonic frequency: 37 kHz Tank capacity: 45 l Tank internal dimensions (W/D/H): 500x300x300 mm for use with S 450 H	POR
FB15150	Fisherbrand* Ultrasonic bath, Elmasonic S 450/(H) Ultrasonic frequency: 37 kHz Tank capacity: 45 l Tank internal dimensions (W/D/H): 500x300x300 mm for use with S 150	POR
FB15081	Fisherbrand* Cover for use with S 10	1344
FB15082	Fisherbrand* Cover for use with S 15	1904
FB15083	Fisherbrand* Cover for use with S 30/ S 40	2128
FB15084	Fisherbrand* Cover for use with S 60	2687
FBEL03283	Fisherbrand* Cover for use with S 70 / S 80	3359
FB15086	Fisherbrand* Cover for use with S100 / S120	3359
FB15087	Fisherbrand* Cover for use with S 180	3583
FB15087	Fisherbrand* Cover for use with S 180	3583
FB15088	Fisherbrand* Cover for use with S300 / S150	4255
FB15093	Fisherbrand* Cover for use with S 450	7615
FB15070	Fisherbrand* Basket for use with S 10	3583
FB15071	Fisherbrand* Basket for use with S 15	3919
FB15072	Fisherbrand* Basket for use with S 30	4591
FB15073	Fisherbrand* Basket for use with S 40	5263
FB15074	Fisherbrand* Basket for use with S 60	6607
FB15075	Fisherbrand* Basket for use with S 70	15341
FB15076	Fisherbrand* Basket for use with S 80	18028
FB15077	Fisherbrand* Basket for use with S 100	13997
FB15078	Fisherbrand* Basket for use with S 120	13997
FB15079	Fisherbrand* Basket for use with S 180	15901
FB15080	Fisherbrand* Basket for use with S 300	18252
FB15094	Fisherbrand* Basket for use with S 450	28554
FB51100	Fisherbrand* Basket for use with S 120	13997
FB15151	Fisherbrand* Basket for use with S 150	18252
FB11040	Fisherbrand* Transparent cover for use with S 15	2464

Prices mentioned are for one unit unless pack size is mentioned in the description

Part No.	Description	Price in Rs.
FB11041	Fisherbrand* Transparent cover for use with S 30 / S 40	3247
FB11042	Fisherbrand* Transparent cover for use with S 60	3695
FB11055	Fisherbrand* Transparent cover for use with S100 / S120	4927
FB11056	Fisherbrand* Transparent cover for use with S300 / S150	7726
FB11035	Fisherbrand* Cover with holes, stainless steel for use with S 10	2687
FB11036	Fisherbrand* Cover with holes, stainless steel for use with S 15	3527
FB11037	Fisherbrand* Cover with holes, stainless steel for use with S30 / S 40	4311
FBEL7093	Fisherbrand* Cover with holes, stainless steel for use with S 60 / S 150	5263
FBEL7038	Fisherbrand* Immersion Basket Hostalen/ diam. = 74mm h =50 mm for use with all models	952
FBEL7054	Fisherbrand* Immersion Basket stainless steel/ diam. =59mm h =60 mm for use with all models	3135
FBEL07545	Fisherbrand* Immersion Basket stainless steel/ diam. =78mm h =60 mm for use with all models	3695
FBEL7039	Fisherbrand* Hostalen tub for use with S 15	1680
FBEL7040	Fisherbrand* Hostalen tub for use with S 30	1904
FBEL7089	Fisherbrand* Acid resistant plastic tub with cover for use with S 30 / S40 S 60	4927
FBEL7090	Fisherbrand* Acid resistant plastic tub with cover for use with S100/ S120	14221
FBEL7091	Fisherbrand* Acid resistant plastic tub with cover for use with S 180	18980
FBEL0080	Fisherbrand* Acid resistant plastic tub with cover for use with S 300	22732
FBEL7078	Fisherbrand* Acid resistant plastic tub with cover for use with S 450	26203
FBEL7036	Fisherbrand* Glass Beaker, diam. =80mm with lid for use with S 10	1277
FBEL7034	Fisherbrand* Glass Beaker, 600ml with lid for use with all models except S 10	1165
FBEL70361	Fisherbrand* Glass Beaker, diam. =80mm without lid for use with S10	952
FBEL70341	Fisherbrand* Glass Beaker, without lid for use with all models except S 10	728
FBEL7076	Fisherbrand* Bottle Holder, stainless steel for use with S 30/ S 40 S 60/ S 70 S 80	2575
FBEL7077	Fisherbrand* Test Tube Holders, stainless steel, 16 Tubes for use with S 30/ S 40 S 60/ S 70 S 80	6719
FBEL07553	Fisherbrand* Clamps for Erlenmeyer flasks 50ml for use with all models	2128
FBEL07555	Fisherbrand* Clamps for Erlenmeyer flasks 250ml for use with all models	2799
FBEL07557	Fisherbrand* Clamps for Erlenmeyer flasks 500ml for use with all models	3359
FBEL07558	Fisherbrand* Clamps for Erlenmeyer flasks 1000ml for use with all models	4031
FBEL07754	Fisherbrand* Cooling coil, stainless steel for use with S 30/ S 60 S 70/ S 80 S 120/ S 180	5487
FBEL09145	Fisherbrand* SRH 4/200 to clean 4 sieves with diam. 200mm at the time for use with S 300	42395
FB11201	Fisherbrand* Sonic bath dual freq.(3 l), Two different ultrasonic frequencies in one unit (frequency can be changed during operation): 37 kHz - for coarse contaminations, for dissolving, blending, dispersing and degassing 80 kHz – noise-reduced operation, ideal for quiet work areas, prolonged ultrasonic cleaning cycle, perfect for the cleaning of hollow spaces, e.g. capillaries for use with P 30	96301
FB11203	Fisherbrand* Sonic bath dual freq. (6 l), Two different ultrasonic frequencies in one unit (frequency can be changed during operation): 37 kHz - for coarse contaminations, for dissolving, blending, dispersing and degassing 80 kHz – noise-reduced operation, ideal for quiet work areas, prolonged ultrasonic cleaning cycle, perfect for the cleaning of hollow spaces, e.g. capillaries for use with P 60	POR
FB11205	Fisherbrand* Sonic bath dual freq.(7 l), Two different ultrasonic frequencies in one unit (frequency can be changed during operation): 37 kHz - for coarse contaminations, for dissolving, blending, dispersing and degassing 80 kHz – noise-reduced operation, ideal for quiet work areas, prolonged ultrasonic cleaning cycle, perfect for the cleaning of hollow spaces, e.g. capillaries for use with P 70	POR
FB11207	Fisherbrand* Sonic bath dual freq.(12 l), Two different ultrasonic frequencies in one unit (frequency can be changed during operation): 37 kHz - for coarse contaminations, for dissolving, blending, dispersing and degassing 80 kHz – noise-reduced operation, ideal for quiet work areas, prolonged ultrasonic cleaning cycle, perfect for the cleaning of hollow spaces, e.g. capillaries for use with P 120	POR
FB11209	Fisherbrand* Sonic bath dual freq.(18 L)Two different ultrasonic frequencies in one unit (frequency can be changed during operation): 37 kHz - for coarse contaminations, for dissolving, blending, dispersing and degassing 80 kHz – noise-reduced operation, ideal for quiet work areas, prolonged ultrasonic cleaning cycle, perfect for the cleaning of hollow spaces, e.g. capillaries for use with P 180	POR
FB11211	Fisherbrand* Sonic bath dual freq.(30 l)Two different ultrasonic frequencies in one unit (frequency can be changed during operation): 37 kHz - for coarse contaminations, for dissolving, blending, dispersing and degassing 80 kHz – noise-reduced operation, ideal for quiet work areas, prolonged ultrasonic cleaning cycle, perfect for the cleaning of hollow spaces, e.g. capillaries for use with P 300	POR

Prices mentioned are for one unit unless pack size is mentioned in the description

Contact: India Toll Free.: 1 800 209 7001 Fax: 022-6680 3001 / 02 e-mail: qfc.customer@thermofisher.com

Part No.	Description	Price in Rs.
FB15083	Fisherbrand* Cover for use with P 30,	2128
FB15084	Fisherbrand* Cover for use with P 60	2687
FBEL3283	Fisherbrand* Cover for use with P 70	3359
FB15086	Fisherbrand* Cover for use with P 120	3359
FB15087	Fisherbrand* Cover for use with P 180	3583
FB15088	Fisherbrand* Cover for use with P 300	3583
FB15072	Fisherbrand* Basket for use with P 30	4591
FB15074	Fisherbrand* Basket for use with P 60	6607
FB15075	Fisherbrand* Basket for use with P 70	15341
FB15078	Fisherbrand* Basket for use with P 120	13997
FB15079	Fisherbrand* Basket for use with P 180	15901
FB15080	Fisherbrand* Basket for use with P 300	18252
FB11037	Fisherbrand* Cover with holes, stainless steel for use with P30 / P 70	4311
FBEL7093	Fisherbrand* Cover with holes, stainless steel for use with P 60	5263
FBEL7038	Fisherbrand* Immersion Basket Hostalen/ diam. = 74mm h =50 mm for use with all models	952
FBEL7054	Fisherbrand* Immersion Basket stainless steel/ diam. =59mm h =60 mm for use with all models	3135
FBEL07545	Fisherbrand* Immersion Basket stainless steel/ diam. =78mm h =60 mm for use with all models	3695
FBEL7040	Fisherbrand* Hostalen tub for use with P 30	1904
FBEL7089	Fisherbrand* Acid resistant plastic tub with cover for use with P 30 / P 60	4927
FBEL7090	Fisherbrand* Acid resistant plastic tub with cover for use with P 120	14221
FBEL7091	Fisherbrand* Acid resistant plastic tub with cover for use with P 180	18980
FBEL0080	Fisherbrand* Acid resistant plastic tub with cover for use with P 300	22732
FBEL7034	Fisherbrand* Glass Beaker, 600ml with lid for use with all models except S 10	1165
FBEL7035	Fisherbrand* Glass Beaker, 100ml with lid for use with P70	2296
FBEL7076	Fisherbrand* Bottle Holder, stainless steel for use with P 30/ P 60 P 70	2575
FBEL7077	Fisherbrand* Test Tube Holderr, stainless steel, 16 Tubes for use with P 30/ P 60 P 70	6719
FBEL07553	Fisherbrand* Clamps for Erlenmeyer flasks 50ml for use with all models	2128
FBEL07555	Fisherbrand* Clamps for Erlenmeyer flasks 250ml for use with all models	2799
FBEL07557	Fisherbrand* Clamps for Erlenmeyer flasks 500ml for use with all models	3359
FBEL07558	Fisherbrand* Clamps for Erlenmeyer flasks 1000ml for use with all models	4031
FBEL07754	Fisherbrand* Cooling coil, stainless steel for use with P30/ P60 P70/ P120 / P180	5487
FBEL09145	Fisherbrand* SRH 4/200 to clean 4 sieves with diam. 200mm at the time for use with P 300	42395
FBEL29450	Fisherbrand* Ultrasonic bath for use with TS 540	89470
FB11089	Fisherbrand* Cover for use with TS 540	4770
FB11085	Fisherbrand* Basket for use with TS 540	8734
FB11301	Fisherbrand* Sonic bath dual freq.(30 Liters),Mains voltage (Vac) 230 - 240 V , Mains connection 1 Ph. / 1N / 1 PE, Power consumption total (W) 2310, Ultrasonic power effective (W) 600, Ultrasonic peak performance max.* (W) 2400, Heating (W) 1700, Unit outer dimensions W/D/H (mm)405 / 449 / 800, Tank internal dimensions W/D/H (mm)330 / 300 / 300, Tank max. volume (L) 30, Tank service volume (L) 22,Basket internal dimensions W/D/H (mm) 255 / 230 / 165,Basket mesh size bottom (mm) 20 / 20 / 2,Basket mesh size site (mm)10 / 10 / 1,6,Basket loading max. (kg)12,Weight (kg)42,Drain,3/4,Material tankstainless steel Material casing stainless steel AISI304 Protection class IP 23 for use with X-tra basic 300	POR
FB11303	Fisherbrand* Sonic bath dual freq.Mains voltage (Vac)400 V , Mains connection 3 Ph. / 1N / 1 PE, Power consumption total (W) 3910, Ultrasonic power effective (W) 1000, Ultrasonic peak performance max.* (W) 4000, Heating (W) 2900, Unit outer dimensions W/D/H (mm) 575 / 479 / 800, Tank internal dimensions W/D/H (mm) 500 / 330 / 350, Tank max. volume (L) 58, Tank service volume (L) 45, for use with X-tra basic 550	POR

Prices mentioned are for one unit unless pack size is mentioned in the description

Part No.	Description	Price in Rs.
FB11305	Fisherbrand* Sonic bath dual freq.Mains voltage (Vac)400 V , Mains connection 3 Ph. / 1N / 1 PE, Power consumption total (W) 6810, Ultrasonic power effective (W) 1000, Ultrasonic peak performance max.* (W) 4000, Heating (W) 2 x 2900, Unit outer dimensions W/D/H (mm) 575 / 479 / 800, Tank internal dimensions W/D/H (mm) 500 / 330 / 500, Tank max. volume (L) 83, Tank service volume (L) 70, for use with X-tra basic 800	POR
FB11307	Fisherbrand* Sonic bath dual freq.Mains voltage (Vac) 400 V , Mains connection 3 Ph. / 1N / 1 PE, Power consumption total (W) 7620, Ultrasonic power effective (W) 1800, Ultrasonic peak performance max.* (W) 7200, Heating (W) 2 x 2900, Unit outer dimensions W/D/H (mm) 676 / 750 / 800, Tank internal dimensions W/D/H (mm) 600 / 600 / 350, Tank max. volume (L) 126, Tank service volume (L) 97, for use with X-tra basic 1200	POR
FB11309	Fisherbrand* Sonic bath dual freq.Mains voltage (Vac) 400 V , Mains connection 3 Ph. / 1N / 1 PE, Power consumption total (W) 9520, Ultrasonic power effective (W) 2000, Ultrasonic peak performance max.* (W) 8000, Heating (W) 3 x 2500, Unit outer dimensions W/D/H (mm) 676 / 750 / 800, Tank internal dimensions W/D/H (mm) 600 / 600 / 450, Tank max. volume (L) 162, Tank service volume (L) 133, for use with X-tra basic 1600	POR
FB11311	Fisherbrand* Sonic bath dual freq.Mains voltage (Vac) 400 V , Mains connection 3 Ph. / 1N / 1 PE, Power consumption total (W) 9520, Ultrasonic power effective (W) 2000, Ultrasonic peak performance max.* (W) 8000, Heating (W) 3 x 2500, Unit outer dimensions W/D/H (mm) 826 / 800 / 800, Tank internal dimensions W/D/H (mm) 750 / 650 / 520, Tank max. volume (L) 255, Tank service volume (L) 215, for use with X-tra basic 2500	POR
FBEL27325	Fisherbrand* Lid Standard for use with 300	5375
FBEL27280	Fisherbrand* Hingend Lid (noise prot.) for use with 300	30234
FBEL27339	Fisherbrand* Basket for use with 300	21836
FBEL26595	Fisherbrand* Lid Standard for use with 550	7279
FBEL26688	Fisherbrand* Hingend Lid (noise prot.) for use with 550	32474
FBEL26584	Fisherbrand* Basket for use with 550	30234
FBEL26595	Fisherbrand* Lid Standard for use with 800	7279
FBEL26688	Fisherbrand* Hingend Lid (noise prot.) for use with 800	32474
FBEL27349	Fisherbrand* Basket for use with 800	36953
FBEL27330	Fisherbrand* Lid Standard for use with 1200	15565
FBEL27297	Fisherbrand* Hingend Lid (noise prot.) for use with 1200	42552
FBEL27350	Fisherbrand* Basket for use with 1200	36953
FBEL27330	Fisherbrand* Lid Standard for use with 1600	15565
FBEL27297	Fisherbrand* Hingend Lid (noise prot.) for use with 1600	42552
FBEL27351	Fisherbrand* Basket for use with 1600	39192
FBEL27332	Fisherbrand* Lid Standard for use with 2500	17357
FBEL27305	Fisherbrand* Hingend Lid (noise prot.) for use with 2500	47031
FBEL27352	Fisherbrand* Basket for use with 2500	47031
FB11160	Fisherbrand* Sonic bath dual freq.TI-H-5 Tank capacity: 4.7 litre Tank internal dimensions (W/D/H): 240 x 130 x 150 mm for use with TI-H 5 MF3	POR
FB11161	Fisherbrand* Sonic bath dual freq.TI-H-5 Tank capacity: 4.7 litre Tank internal dimensions (W/D/H): 240 x 130 x 150 mm for use with TI-H 5 MF2	POR
FB11163	Fisherbrand* Sonic bath dual freq.TI-H-10 Tank capacity: 8.6 litre Tank internal dimensions (W/D/H): 300 x 240 x 150 mm for use with TI-H 10 MF3	POR
FB11164	Fisherbrand* Sonic bath dual freq.TI-H-10 Tank capacity: 8.6 litre Tank internal dimensions (W/D/H): 300 x 240 x 150 mm for use with TI-H 10 MF2	POR
FB11166	Fisherbrand* Sonic bath dual freq.TI-H-15 Tank capacity: 12.2 litre Tank internal dimensions (W/D/H): 300 x 240 x 200 mm for use with TI-H 15 MF3	POR
FB11167	Fisherbrand* Sonic bath dual freq.TI-H-15 Tank capacity: 12.2 litre Tank internal dimensions (W/D/H): 300 x 240 x 200 mm for use with TI-H 15 MF2	POR
FB11169	Fisherbrand* Sonic bath dual freq.TI-H-20 Tank capacity: 16,8 litre Tank internal dimensions (W/D/H): 330 x 300 x 200 mm for use with TI-H 20 MF3	POR
FB11170	Fisherbrand* Sonic bath dual freq.TI-H-20 Tank capacity: 16,8 litre Tank internal dimensions (W/D/H): 330 x 300 x 200 mm for use with TI-H 20 MF2	POR
FB11301	Fisherbrand* Sonic bath dual freq.FB11301 - Bath X-tra basic 300 ultrasonic, floor standing multi-frequency 25/45kHz, stainless steel 30L tank volume 405mm x 449mm x 800mm (w x d x h) heater power 1 x 1700W Fisherbrand for use with TI-H 25 MF2	POR
FB11305	Fisherbrand* Sonic bath dual freq.FB11305 - Bath X-tra basic 800 ultrasonic, floor standing multi-frequency 25/45kHz, stainless steel 83L tank volume 575mm x 479mm x 800mm (w x d x h) heater power 2 x 2900W Fisherbrand for use with TI-H 80 MF2	POR
FB11307	Fisherbrand* Sonic bath dual freq.Bath X-tra basic 1200 ultrasonic, floor standing multi-frequency 25/45kHz, stainless steel 126L tank volume 676mm x 750mm x 800mm (w x d x h) heater power 2 x 2900W Fisherbrand for use with TI-H 115 MF2	POR
FB11309	Fisherbrand* Sonic bath dual freq.Bath X-tra basic 1600 ultrasonic, floor standing multi-frequency 25/45kHz, stainless steel 162L tank volume 676mm x 750mm x 800mm (w x d x h) heater power 3 x 2500W Fisherbrand for use with TI-H 160 MF2	POR

Prices mentioned are for one unit unless pack size is mentioned in the description